



STANDING FRAME

**GENIUS<sup>®</sup>**



USER MANUAL





Manufacturer takes responsibility for proper usage of the device only when the product was purchased from qualified **Akces-MED Ltd.** employee or at specialized medical store!



**Akces-MED Ltd.** reserve the right to introduce technical and trading changes in the content of the instruction without warning.



Product marked with the CE sign is Class I medical device, of rule 1, consistent with the requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.



**Akces-MED Ltd.** Company introduced quality management system ISO 13485 regarding production, sales and service. The system is certified by Det Norske Veritas.





## CONTENTS

1. The purpose of user manual.....	5
2. Indications and contrindications for use.....	5
3. General characteristic.....	7
3.1 Standard equipment.....	7
4. Size table.....	8
5. Description how to use standing frame.....	9
5.1 Chest and pelvic belt adjustment.....	10
5.1.1 Height adjustment of the belts.....	10
5.1.2 Setting the position and adjusting the width of the chest and pelvic belts.....	11
5.1.3 Depth and rotation adjustment of the belts.....	14
5.1.4 Belt support.....	16
5.1.5 Chest and pelvic belts.....	17
5.2 Tray.....	18
5.2.1 Tray height adjustment.....	18
5.2.2 Tray depth adjustment.....	19
5.2.3 Tray tilt angle adjustment.....	20
5.3 Knee supports.....	21
5.3.1 Knee supports height adjustment.....	21
5.3.2 Knee supports spacing adjustment.....	22
5.3.3 Front-back knee supports adjustment.....	23
5.3.4 Knee supports angle adjustment.....	24
5.3.5 Front-back knee supports angle adjustment.....	25
5.4 Sandals with the adduction and abduction function, equipped with feet stabilizing belts.....	26
5.4.1 Sandals spacing adjustment.....	26
6. The use of standing frame.....	28
6.1 The preparation of user.....	28
6.2 The preparation of standing frame.....	29
6.3 Upright position.....	30
7. Safety rules.....	32
8. Terms of warranty.....	34
9. Maintenance and cleaning.....	35
10. Labelling.....	39
10.1 Symbols.....	39
10.2 Sample data plate.....	39





## INTRODUCTION

The **GENIUS®** standing frame has been created for the most demanding customers. Its well thought-out design gives you infinite adjustment possibilities of all the independent elements in every plane so the standing frame can be perfectly adjusted depending on the user's needs.

Every element is thoroughly designed and the materials used are of very high quality – all these factors make the standing frame comfortable and safe for the user. Rehabilitation of people with spinal, pelvic or knee joint dysfunctions does not have to be complicated – with the **GENIUS®** standing frame it is much more effective and its effects are long-lasting.

*Remember to read this user manual carefully. You will find not only useful information but also guidelines regarding safe and effective use of the product.*

*If you have any questions or remarks, do not hesitate to contact with us:*

**Medical consultant:** mobile +48 516 183 883

**Akces-MED. Sp. z o.o.:** tel. +48 17 864 04 77

**e-mail:** [export@akces-med.com](mailto:export@akces-med.com)

***Thank you for purchasing our product!***



Before using the stander, read this user manual which is essential for proper and safe usage of the stander. Remember that complying with guidelines is highly important. It assures safety and extends durability as well as aesthetics of the device.







## 1. THE PURPOSE OF USER MANUAL



**The present user manual contains basic information which is essential for :**

- preparation before use,
  - proper usage,
- maintenance, cleaning and guarantee conditions.

User manual should be kept in easily accessible place. It is designed for people who look after the disabled , and for both doctors as well as physiotherapists who operate the device. Before using the **GENIUS®** standing frame, user is obliged to familiarize himself/herself with user manual.

## 2. INDICATIONS AND CONTRINDICATIONS FOR USE

The standing frames are intended for people with orthopedic, neurological and cardiological diseases and for people after surgery who are unable to adopt and maintain a standing position on their own as a result of prolonged immobilization. The device is an irreplaceable aid in the rehabilitation process, both in terms of movement and social integration. It is great help in everyday activities such as rehabilitation but also during playing or studying. The selection of the appropriate size of the device depends on medical indications and anatomical parameters, regardless of the age of the user.

**It is recommended for users in the following cases:**

- cerebral palsy – various types and forms,
- brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology,
- demyelinating disease (e.g. multiple sclerosis – SM),
- meningocele,
- muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis),
- genetic disorders,
- egeneration and neurological diseases,





- other condition connected with muscles weakness, paralysis and
- musculoskeletal disorders,
- after head and spinal cord injuries and spine trauma with paralysis and paresis.

## **Contraindications:**

While standing upright, the user who is in a standing position should be closely monitored. The therapy should be discontinued in case of such side effects as: turning pale, sweating or fainting, joint strain and the resulting deformities of the body. Further proceeding should be consulted with the medical specialist or a physiotherapist who prescribed the treatment. The critical condition of the user is an absolute contraindication to standing upright.

Relative contraindications for upright standing are: arterial hypertension at rest above 200/110 mmHg, pulmonary hypertension, cardiac arrhythmias, unstable asthma, unstable metabolic disorders, moderate heart valve disease.

- Any condition, such as a healing fracture or severe osteoporosis, that makes it impossible to bear the user's own bodyweight.
- Bone fracture with incomplete union, unless approved by a doctor.
- Postoperative restrictions on joint loading.
- Pain while standing in the standing frame (caused by subluxation or dislocation of the hip joint, overstretching of soft tissues or excessive pressure). Subluxation or dislocation of the hip joint in the absence of pain is not a contraindication.

People with significant flexion in the hip or knee joints, as well as with contractures in the ankle joint, may not be able to use the standing frame. Make sure that the position of the user in the standing frame does not cause excessive pressure or stretching of affected areas of the body. In such cases, standing upright may be contraindicated until a greater range of motion is achieved in the contracted joints.

In case of people with impaired cardiovascular system (arterial hypertension at rest above 200/110 mmHg, pulmonary hypertension, cardiac arrhythmias, moderate heart valve disease); respiratory system (unstable asthma) and unstable metabolic disorders, it is necessary to monitor the user while standing in order to prevent serious complications.





## 3. GENERAL CHARACTERISTIC



### 3.1 STANDARD EQUIPMENT

- ① Frame
- ② Mobile base with wheels with brakes
- ③ Sandals with the adduction and abduction function, equipped with feet stabilizing belts
- ④ Independently adjustable knee supports
- ⑤ Pelvic belt with pelvic rotation function
- ⑥ Chest belt with chest rotation function
- ⑦ Tray





## 4. SIZE TABLE



SIZE TABLE

Symbol	Measurement	Unit	Size 1	Size 2	Size 3
	indicative user height*	[cm]	75 - 100	95 - 125	120 - 150
A/A <sub>1</sub>	chest/pelvic support width	[cm]	10 - 24	12 - 24	20 - 33
C	chest support height**	[cm]	55 - 67	72 - 89	83 - 107
C <sub>1</sub>	pelvic support height**	[cm]	45 - 55	56 - 69	67 - 87
B/B <sub>1</sub>	chest/pelvic support perimeter	[cm]	43 - 67	50 - 77	66 - 100
C <sub>2</sub>	tray height**	[cm]	54 - 70	64 - 90	73 - 100
G	knee supports width	[cm]	6 - 25	6 - 30	6 - 35
G <sub>1</sub>	knee supports height**	[cm]	25 - 35	23 - 43	37 - 53
K	footer length	[cm]	17	20	23
	max. user weight	[kg]	35	45	55

Weight and dimensions		Size 1	Size 2	Size 3
Width	[cm]	56	60	66
Length	[cm]	72	80	88
Height	[cm]	66	83	95
Weight	[kg]	20	24	27

**Note:** Measurement tolerance range is +/- 1 cm / ½ in.

\* Height cannot be the main parameter while choosing the product size.

\*\* Height measured from the foot adjustment.





## 5. DESCRIPTION HOW TO USE STANDING FRAME

### KEY



**Info icon** - appears in the manual when you need to read the additional information.



**Info icon** - appears in the manual when you need to use the tools indicated in the assembly description.



**Warning icon** - appears in the manual when you need to take special precautions.



## 5.1 CHEST AND PELVIC BELT ADJUSTMENT

### 5.1.1 HEIGHT ADJUSTMENT OF THE BELTS



**i**

Use knobs **A** to adjust the height of the chest belt (**fig. 2 - 3**). Before adjusting the height, loosen both knobs.

**1**



**2**




**3**





Use knobs **(B)** to adjust the height of the pelvic belt (**fig. 4 - 5**). Before adjusting the height, loosen both knobs.



## 5.1.2 SETTING THE POSITION AND ADJUSTING THE WIDTH OF THE CHEST AND PELVIC BELTS





**i**

Use knobs **C** to adjust the width of the chest belt (**fig. 10 - 11**)  
and to position the belt laterally (**fig. 12 - 13**).

Use knobs **D** to adjust the width of the pelvic belt (**fig. 10 - 11**)  
and to position the belts laterally (**fig. 12 - 13**).







## WIDTH ADJUSTMENT OF THE BELTS



## LATERAL ADJUSTMENT OF THE BELT





## 5.1.3 DEPTH AND ROTATION ADJUSTMENT OF THE BELTS

i

Use knobs **C** to adjust the depth (fig. 16 - 17) and rotation (fig. 21 - 24) of the chest belt.

Use knobs **D** to adjust the depth (fig. 16 - 17) and rotation (fig. 21 - 24) of the pelvic belt.



### DEPTH ADJUSTMENT OF THE BELT





18



19



20

## ROTATION ADJUSTMENT OF THE BELT



21



22





## 5.1.4 BELT SUPPORT



Both the chest and pelvic belts have a special support. The support helps to adjust the belts to the user's body.





## 5.1.5 CHEST AND PELVIC BELTS



The chest and pelvic belts can be adjusted to the circumference of the chest and pelvic with the use of the Velcro fastening. The fastening is located at the front and back of the belts.





## 5.2 TRAY

### 5.2.1 TRAY HEIGHT ADJUSTMENT



**31**

**i**

Use knobs **E** to adjust the height of the tray (**fig. 31 - 33**). Before adjusting the height, loosen both knobs.



**32**



**33**





## 5.2.2 TRAY DEPTH ADJUSTMENT



34



Use knobs **F** to adjust the depth of the tray (**fig. 34 - 36**). Before adjusting the depth, loosen both knobs.



35



36





## 5.2.3 TRAY TILT ANGLE ADJUSTMENT



Use knobs **G** to adjust the tilt angle of the tray (**fig. 37 - 40**).  
Before adjusting the angle, loosen both knobs.







## 5.3 KNEE SUPPORTS

### 5.3.1 KNEE SUPPORTS HEIGHT ADJUSTMENT



i

Use knobs **(K)** to adjust the height of the knee supports (**fig. 41 - 43**). Each knee support can be adjusted independently.






## 5.3.2 KNEE SUPPORTS SPACING ADJUSTMENT



Use knobs **(H)** to adjust the spacing of the knee supports (**fig. 44 - 48**).

You can adjust the knee supports according to the deformations within the knee joints – valgus or varus. Each knee support can be adjusted independently.





## 5.3.3 FRONT-BACK KNEE SUPPORTS ADJUSTMENT



i

Use knobs **J** to adjust the knee supports in the front-back plane (fig. 49 - 53).

You can adjust the knee supports depending on the needs and therapeutic indications, e.g. existing flexion contractures. Each knee support can be adjusted independently.





## 5.3.4 KNEE SUPPORTS ANGLE ADJUSTMENT



54



Use knobs **I** to adjust the rotation of the knee supports (fig. 54 - 58).

You can set the knee supports at the proper angle to achieve the desired positioning of the knee joint area.


Each knee support can be adjusted independently.




55



56



57



58





### 5.3.5 FRONT-BACK KNEE SUPPORTS ANGLE ADJUSTMENT



With the additional joint, the knee supports automatically adjust to the user.





## 5.4 SANDALS WITH FEET STABILIZING BELTS



The user of the stander must always wear proper footwear which provides stabilization of the ankle joints to minimize the risk of injury.




### 5.4.1 SANDALS SPACING ADJUSTMENT



Use knobs  to adjust the spacing of the sandals (**zdj. 64 - 66**).

Sandals can be set in a selected position conducive to the correct upright standing







## 6. THE USE OF STANDER



Before using the **GENIUS®** stander, read this user manual.  
Before each use, prepare the user and the device!

### 6.1 PREPARATION OF THE USER

**Before using the product, it is essential to make sure that:**

- user's outfit is light and comfortable sporty clothes are advisable,
- shoes should be comfortable and, moreover, they should stabilize ankle,
- due to hygienic reason and to avoid abrasion it is advisable for the user to always wear clothes when using the device, skin must not be naked especially in places of regular contact with upholstery,
- drains, bags, catheters etc. if applied, have to be secured to function properly and not to disturb during using the product,
- if it is necessary for a patient to wear orthopedic supply (collar, corset, orthosis, shoes) if it is recommended by doctor.



Outfit should protect user's skin from direct contact with upholstery as well as other parts of equipment.

Only doctor can make a decision concerning using any kind of orthopedic supply!







## 6.2 THE PREPARATION OF STANDING FRAME

### 1. How to prepare the standing frame for use:

- assemble the standing frame;
- place the standing frame on an even and flat surface, away from sources of fire and high temperature;
- make sure that all the elements are undamaged, properly mounted and stable;
- push the brakes;
- adjust the standing frame to the user's size and needs;
- unfasten the chest and pelvic belts in order to place the user in the standing frame;
- fasten the elements which stabilize the user - the chest and pelvic belts;
- mount the tray.



**Before using the standing frame, make sure that all the belts and other elements have been properly adjusted, and all the knobs and fastenings of the adjustable elements are firmly tightened.**

**The standing frame will fulfil its function properly only when it is adjusted to the size and needs of the user. The most important thing, while adjusting the standing frame to the user, is to provide the user with the maximum feeling of safety.**

**The improper use of the standing frame may lead to strengthening the incorrect posture patterns.**





## 6.3 UPRIGHT POSITION

**It is necessary to:**

1. Prepare a place for upright standing. It is important to choose a room without any dangers (such as: strewn toys or other objects).
2. Move a standing frame closer to a place in which the user is situated.
3. Block castors from a standing frame, wheel chair or other device in which the user is sitting.
4. Move the user and place him in a standing frame. We ought to make sure that all alignments are preserved. Implement corrections if necessary.
5. Fasten belts stabilizing chest, hips and foot.
6. Make sure that all belts and precautions were correctly adjusted and fastened.



**The length of time for upright position process is different for each user. It is recommended to extend this time during the initial period (adaption period for the change of position).**

**Either doctor or physiotherapist can make decision concerning the length of upright position process. Upright position process has to be always performed under the supervision of trained person!**

**It is recommended to divert user's attention during upright position process by:**

- manipulation games,
- educational games,
- touch stimulation and others.



**It is highly important to adjust the device to the user before using it!**

**Make sure that all belts as well as fastenings from adjustable elements are correctly adjusted and fastened!**

**It is forbidden to leave the patient unattended!**





**Incorrect adjustment of the device can lead to side effects such as: fainting, body deformities, overloading joints. It is highly advisable to adjust the device before every use of a standing frame. This adjustment has to be done only by a trained person or therapist. If any of these side effects appear, it is highly important to stop the upright position process, take out the user from the device and lie him/her down.**

**If these undesirable side-effects still occur it is necessary to contact doctor.**

**If these side effects disappeared, remember to inform doctor so as to make decision concerning further upright position process.**

### **The ending of upright position:**

1. It is essential to prepare a place in which the user will be taken out from a standing place. It is important to choose a room without any dangers.
2. Move a standing frame closer to a place in which the user will be moved.
3. Block castors from a standing frame, orthopaedic bed or other device in which the user will be situated.
4. Unfasten belts which stabilize foot, hips and chest.
5. Make sure that all belts and precautions do not block the user.
6. Move the user and place him on a prepared device.





## 7. SAFETY RULES



Stander **GENIUS®** has to be used in compliance with guidelines in this user manual.  
Before each use, prepare the user and the device.

It is highly important to consult either doctor or physiotherapist in order to make right decision concerning the choice of **GENIUS®** and its size. It is also essential to familiarize with the guidelines described in user manual before using the device.

1. Before use it is necessary to:

- make sure that there are no damages or different faults which could endanger user's life; check if all parts are working and are correctly mounted and, moreover, check if they are not broken;
- check if all screws, knobs and all glued, sewed or tapped parts are tightened;
- check the condition of the upholstery for any signs of damage;
- check if brakes are operating correctly;
- check if buckles are operating correctly;
- always use safety belts, Velcro tapes should overlap min. 8cm;
- adjust the device according to user's measurements and needs e.g. the adjustment of pelottes etc.;
- make sure that all movable parts and knobs for adjustments are stable and locked;
- secure the user with available stabilizing aids of the standing frame. It is highly important to secure patient because it safeguards him/her from falling out.

2. It is forbidden to use standing frame for different purposes that it is intended





3. It is forbidden to use standing frame when it is damaged or any part is missing.
4. It is forbidden to leave the user unattended.
5. Stroller should be used only when there is someone who is trained and knowledgeable.
6. It is forbidden to leave standing frame and the user on unstable ground unattended, even if the brake is engaged.
7. It is forbidden to either go up or down the stairs; it doesn't matter if the user is placed in the device or not.
8. It is important to engage brake during getting into and getting out the standing frame.
9. Make sure that all limbs are not close to spokes or other adjustable elements.
10. The standing frame should be stored and used away from sources of fire or high temperature. Moreover, standing frame shouldn't be stored in rooms with high humidity
11. Secure standing frame during transportation in order to avoid any possible damages. Manufacturer doesn't take responsibility for damages which appeared due to improper preparation of the device during transportation.
12. Before using the standing frame, remove all paper and plastic packaging. The packaging is not a toy! Keep away from children and animals.



**Incorrect product adjustment can lead to such side effects as fainting, body deformities, joint overload, etc. If the abovementioned symptoms occur, stop using the product immediately, take the user out of the product and lay them down. If the undesirable symptoms persist, contact a physician immediately. If the symptoms have disappeared, inform the attending physician about this fact in order to make a decision about the further rehabilitation process.**

**Serious incidents involving this device should be reported both to manufacturer and to the competent authority in your country.**





**The device contains small parts which children may swallow.**

**Do not use a product that shows any signs of damage, destruction or wear.**

**Only a fully functional product, used in accordance with the instructions for use, is safe for the user. If the product is damaged, discontinue use and contact us with the manufacturer.**

**It is forbidden to make any changes concerning device's structure due to the fact that it is dangerous and it results in losing warranty.**

**The weight of user should not be bigger than maximal weight described on page 8 where you will find information regarding this issue.**

**This device can be only used inside the buildings!**

## 8. TERMS OF WARRANTY

- **Akces-MED Ltd.** sells products free from manufacturing defects or faults.
- Manufacturer gives the purchaser 24-month guarantee from the date of purchase. All defects revealed during warranty period will be repaired free of charge.
- All defects which appear after warranty period can be repaired only when the user will pay for it. This rule doesn't apply to refunded products - user can ask National Health Refund for refund.
- The stroller is covered by a 36-month warranty on the territory of Canada – the details of the rights resulting from the extended warranty are regulated by the warranty card.
- Manufacturer claims that the product can be used for 48-months. User can use this product after this period only when he/she will send the device to manufacturer so as to service it. Next, manufacturer determines again how long the product can be used and, moreover, the date of the next service will be determined.
- If product is damaged, please contact our dealer or **Akces-MED Ltd.** customer service (ph. +48 17 864 04 77, e-mail: [export@akces-med.com](mailto:export@akces-med.com)).





## 9. MAINTENANCE AND CLEANING

The standing frame **GENIUS®**, which is a medical device, should always be kept clean and should be used in accordance with the manufacturer's instructions. The standing frame should be kept in a clean, dry and ventilated room. Maintenance and cleaning must be done at least once a month. Periodic cleaning and maintenance of the device is essential to keep its proper functioning and safety.

### **Recommendations for cleaning the upholstery and other elements of the device:**

1. Wash with warm water using a soft cloth and mild cleaners.
2. Make sure the cloth is not too wet, as this may enlarge the stain and damage the upholstery or foam inside the upholstery.
3. Remember to remove excess moisture from the upholstery before the next cleaning.
4. The upholstery must not get wet during cleaning.
5. You must not soak the upholstery.!
6. After cleaning the stains, the upholstery should be rinsed with lukewarm water and dried with a clean cloth.
7. Do not use hard brushes for cleaning the upholstery.
8. Do not clean the upholstery with strong chemicals, e.g. solvents or chlorine-based caustics.
9. Do not clean, replace any elements or maintain the standing frame during use.
10. Metal and plastic elements should be cleaned by hand with suitable cleaners, after having read the recommendations of the manufacturer of the given cleaner
11. It is possible to disinfect metal, plastic, wooden and upholstery elements of the product with a generally available disinfectant. In order to disinfect the selected element, spray its surface with a disinfectant generally available on the market and wipe with a dry, clean cloth. After disinfection, wait until the element is completely dry, especially in the case of upholstery. Before using the disinfectant, first check in an invisible place whether the product will not discolor or damage the disinfected surface.





Remember that upholstery has to be always dry before using!

**The manufacturer of the standing frame, Akces-MED Sp. z o. o., is not responsible for any damages caused by using unsuitable cleaners.**

### Before using it is necessary to:

- clean wheels;
- check if base as well as wheels are correctly mounted; check if brake stops turning wheel;
- check if standing frame construction is free of breaks or any kinds of deformities;
- check if all movable parts are working and all screws, nuts, knobs etc. are correctly tightened;
- check if belts are not frayed; check if buckles and clamps work correctly.



**Manufacturer takes responsibility only for hidden faults or damages which appeared despite the fact that user complied with the guidelines.**

In case of any damages, it is essential to contact the person who sold you the product.







## TECHNICAL INSPECTION

Technical inspections should be performed regularly. Before servicing, the product should be thoroughly cleaned to reveal any hidden damage.

Frequency	What should be checked? What should be done?		Remarks
Before each use.	General condition of the product.	<b>Check:</b> <ul style="list-style-type: none"> <li>the entire product for any visible damage, cracks or mechanical defects,,</li> <li>if the product complete,</li> <li>if all bolts, screws, knobs, nuts and other elements of the construction are complete and tightened..</li> </ul>	In case of any damage or defects, discontinue the use of the product and contact your Supplier/Distributor or the Manufacturer.
	Wheels and brakes.	<b>Check:</b> <ul style="list-style-type: none"> <li>if the wheels are properly mounted and if they are not damaged,</li> <li>if the brake is blocking the spinning wheel,</li> <li>if the wheels are not dirty. If necessary, clean the wheels.</li> </ul>	In case of any damage or defects, discontinue the use of the product and contact your Supplier/Distributor or the Manufacturer.
	Fastening of the construction elements.	<b>Check all the mechanisms of the product and its moving parts. Check in particular:</b> <ul style="list-style-type: none"> <li>if the seat is properly mounted (if applicable),</li> <li>if the wheels are properly mounted (if applicable),</li> <li>if the elements securing the user are properly mounted and if they work well (if applicable),</li> <li>if all the moving parts are complete, work properly, do not jam,</li> <li>if the adjustment mechanisms of the following elements work properly (if applicable): headrest, backrest, seat, abduction block, armrests, footrests.</li> </ul>	In case of any damage or defects, discontinue the use of the product and contact your Supplier/Distributor or the Manufacturer.
	Upholstery elements.	<b>Check:</b> <ul style="list-style-type: none"> <li>if there is any damage, tears, rips, loose seams or deformations on the upholstery elements,</li> <li>if the stabilizing straps, vests or safety belts are not frayed or twisted and if the buckles work properly,</li> <li>if the stabilizing straps, vests or safety belts work correctly and do not unfasten under pressure or when pulled.</li> </ul>	In case of any damage or defects, discontinue the use of the product and contact your Supplier/Distributor or the Manufacturer.
Once a week.	Product cleaning.	If the frame and upholstery get dirty (dusty, muddy, etc.) during everyday use, wipe them with a damp cloth without detergents. Remove dirt (e.g. hair, food remains, etc.) from the moving parts of the product.	The product should be cleaned when necessary, but at least once a week. Do not use any chlorine based cleaners or methylated spirit. Before using the product, check if the upholstery is dry.





Frequency	What should be checked? What should be done?		Remarks
Once a month.	Wheel axles, quick connectors (if applicable).	<b>Check:</b> <ul style="list-style-type: none"> <li>if the wheels are rolling and if they are not loose. To maintain the wheels, apply a small amount of a preservative-lubricant between the axle and the wheel hub. Check if you can easily remove and attach the wheels.</li> </ul>	Clean the wheels when dirty. In case of any damage or defects, discontinue the use of the product and contact your Supplier/Distributor or the Manufacturer.
	Product construction connectors (welding, threading, soldering)	<b>Check:</b> <ul style="list-style-type: none"> <li>the condition of the product frame in the places of joining (welds, threads, soldering).</li> </ul>	In case of any damage or defects, discontinue the use of the product and contact your Supplier/Distributor or the Manufacturer.
	Tires/innertubes (if applicable).	<b>Check:</b> <ul style="list-style-type: none"> <li>tire wear condition. Check for visible damage. Check tire pressure (if applicable).</li> </ul>	If the tire pressure is too low, inflate them to the correct level. In case of any damage or defects, discontinue the use of the product and contact your Supplier/Distributor or the Manufacturer.
Once every 6 months.	Construction moving parts.	<b>Check:</b> <ul style="list-style-type: none"> <li>if the moving parts of the construction have not been damaged. Maintain them using a proper cleaning agent.</li> </ul>	In case of any damage or defects, discontinue the use of the product and contact your Supplier/Distributor or the Manufacturer.
	Product labeling.	<b>Check:</b> <ul style="list-style-type: none"> <li>if the marking and labels are readable.</li> </ul>	In case of any damage or defects, discontinue the use of the product and contact your Supplier/Distributor or the Manufacturer.
After the life time of the product finishes.	Detailed product check.	Check the product.	The manufacturer specifies the life time of the product, which is given in the user manual. After the given period, the user may use the product only if it is sent back to the manufacturer to be checked first. After that, the manufacturer will re-define the life time of the product and the date of the next check.

The list of all spare parts is available for download after selecting a product on [www.akces-med.com](http://www.akces-med.com).





## 10. LABELLING

### 10.1 SYMBOLS

Declaration that the product meets the essential requirements

Catalog number

Serial number

Please read the manual instructions before use

You need to take special precautions

Manufacturer

Production date

The product does not contain latex

Store in a dry place

Careful! Delicate content

Protect from sunlight

Storage temperature limitation

Relative humidity storage

Dry vertically

Do not chlorinate

Do not iron

Do not dry in a drum dryer

Do not bleach

The product is intended for indoor use

Prohibition of use without fastened safety belts

Prohibition to leave the patient without care during use of the product

You need to read the additional information

You need to use the tools indicated in the assembly description

### 10.2 SAMPLE DATA PLATE

**QR CODE**  
after scanning, it directs you to the product's website

**MANUFACTURER ADDRESS**  
AKCES-MED Sp. z o.o.  
Jasionka 955B  
36-002 JASIONKA  
POLAND

**PRODUCT NAME PL**  
GENIUS 3  
PIONIZATOR STATYCZNY

**PRODUCT NAME EN**  
GENIUS 3  
STANDING FRAME

**PRODUCT BARCODE**  
5 906395 195862

**MAXIMUM PERMISSIBLE LOAD ON THE PRODUCT**  
MAX LOAD: 55 kg

**DATE OF PRODUCTION**  
07.2021

**SERIAL NUMBER IN BARCODE FORM**  
SN GNS 033/U/3/BA

**SERIAL NUMBER**

**LABELLING look up**  
MD, Manual icon, CE




## MANUFACTURER



**Akces-MED Sp. z o.o.**

Jasionka 955B  
36-002 Jasionka  
tel. +48 17 864 04 70  
e-mail: [export@akces-med.com](mailto:export@akces-med.com)  
BDO: 000031619



Authorised Dealer



**akces-med.com**

Date of issue: 08.01.2020  
The last update: 27.07.2021  
©Sławomir Wroński 2021