



SPECIAL NEEDS STROLLER  
NOVA EVO<sup>®</sup>



USER MANUAL





Manufacturer takes responsibility for proper usage of the device only when the product was purchased from qualified Akces-Med employee or at specialized medical store!

AKCES-MED. Ltd. reserve the right to introduce technical and trading changes in the content of the instruction without warning.

**MEDICAL  
DEVICE**  
CLASS 1





**WYRÓB  
MEDYCZNY**  
KLASA 1

Medical product first class according to the Regulation 93/42/EWG of European Parliament and Council regarding medical products. For this product the producer issued relevant declaration of conformity.



Akces-Med. Sp. z o.o. Company introduced quality management system ISO 13485 regarding production, sales and service. The system is certified by Det Norske Veritas.





## CONTENTS

1. The purpose of user manual.....	6
2. Indication for use.....	6
3. General characteristics of stroller.....	7
3.1. Standard equipment.....	7
3.2. Accessories.....	8
4. Size table.....	10
5. How to use the stroller.....	11
5.1 Floding the stroller.....	12
5.2 Unfloding the stroller.....	13
5.3 Footrest.....	15
5.3.1 Footrest tilt angle adjustment.....	15
5.3.2 Footrest length adjustment.....	15
5.3.3 Footplate tilt angle adjustment.....	16
5.4 Wheels.....	16
5.4.1 Removing rear inflatable wheels.....	16
5.4.2 Removing front swivel PU wheels.....	17
5.4.3 Mounting front inflatable wheels / smartraction system.....	18
5.4 Central brake.....	20
5.5 Adjustable backrest.....	21
5.6 Tilt-in-space function.....	22
5.7 Headrest.....	23
5.8 5-point safety belts.....	24
5.9 Upholstery.....	28
5.9.1 Flex upholstery.....	28
5.9.2 Waterproof upholstery HYDRO.....	31
5.10 Safety rail.....	33
6. Accessories.....	35
6.1 NVE_627 NOVA module.....	35
6.2 NVE_625 HOME module.....	36
6.3 NVE_128 Abduction block.....	38





6.4 NVE_405 Flooding canopy.....	38
6.5 NVE_445 Footmuff.....	39
6.6 NVE_114 H Harness / NVE_130 4 point safety vest.....	40
7. The use of stroller.....	45
7.1 Preparation of the patient.....	45
7.2 Preparation of the stroller.....	46
7.3 How to use vest.....	47
8. Safety rules.....	50
9. Term of warranty.....	52
10. Maintenance and cleaning.....	53
11. Labelling.....	55
10.1 Symbols.....	55





## INTRODUCTION

Special needs stroller **NOVA EVO®** is a great combination of a stroller which is perfect both during home rehabilitation and long walks outside. It is especially recommended for spastic users as well as those requiring additional stabilization. The deep seat provides a comfortable position in the stroller, which helps to correct posture defects and significantly prevents their deepening.

Numerous adjustments and a wide range of the accessories make it possible to adapt the stroller to the individual needs of the user. An adjustable footrest and the 3-step tilt-in-space function of the seat guarantee a comfortable and convenient position in the stroller.

**NOVA EVO®** is equipped with front swivel PU wheels and rear inflatable wheels with central brake. The stroller is easy to fold thanks to its construction and the compact dimensions make it fit even into the smallest car trunks.

If you have any questions or remarks, do not hesitate to contact with us:

**Medical consultant: mobile +48 508 382 509**

**Akces-MED Sp. z o.o.: ph. +48 17 864 04 77**

**e-mail: [export@akces-med.com](mailto:export@akces-med.com)**

***Thank you for purchasing our product***



Before using the product, user is obliged to familiarize himself/herself with user manual. Remember that complying with guidelines is highly important. It assures safety and extends durability as well as aesthetics of the device.





## 1. THE PURPOSE OF USER MANUAL



**The instruction manual contains basic information necessary to properly:**

- prepare the stroller before use
- use the stroller
- maintain and clean the stroller

This user manual should be kept in an easily accessible place. The manual is intended for people who look after disabled users, as well as for doctors and physiotherapists handling the stroller. Before using the stroller for the first time, the user (or caregiver) is obliged to familiarize himself/herself with the user manual

## 2. INDICATIONS FOR USE

The device is recommended for people in the following cases:

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis - SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.





## 3. GENERAL CHARACTERISTICS OF STROLLER

### 3.1 Standard equipment

- ① Foldable frame
- ② Removable rear inflatable wheels
- ③ Front swivel PU wheels
- ④ Adjustable rear wheelbase SmarTrack
- ⑤ Adjustable footrest with 2-point feet stabilizing belts
- ⑥ Seat cushion
- ⑦ Push handle
- ⑧ 5-point safety belts
- ⑨ Adjustable headrest
- ⑩ Safety rail
- ⑪ Central brake
- ⑫ Flex upholstery
- ⑬ Waterproof HYDRO™ upholstery
- ⑭ Tilt-in-space adjustment





### 3.2 Accessories



NVE\_001

LED lights



NVE\_002

Belt for moving patients AM-P



NVE\_101

Thigh abduction belts



NVE\_107

Pelvic belt



NVE\_114

H Harness



NVE\_126

Trunk belt



NVE\_128

Adjustable abduction block



NVE\_130

4 points safety vest



NVE\_134

Narrowing pads 6 cm



NVE\_137

Narrowing pads 10 cm



NVE\_147

5 points feet stabilizing belts



NVE\_161

Foot and ankle stabilizing sandals



NVE\_162

Dynamic ankle stabilizers (FP-22)



NVE\_403


Tray



NVE\_405


Floding canopy







## 4. SIZE TABLE



### SIZE TABLE

Symbol	Measurement	Unit	Size 1	Size 2	Size 3
W	Indicative user height*	[cm]	100-130	120-160	140-180
A	Seat width	[cm]	31	36	41
L	Seat depth	[cm]	32	38	44
C1	Backrest height	[cm]	62	70	84
C	Headrest height	[cm]	43-55	54-66	64-79
E	Backrest width	[cm]	31	36	41
F	Footrests length	[cm]	25-43	30-50	32-56
	Footer length	[cm]	21	24	27
	Push handle height	[cm]	96	96	96
X	Backrest tilt angle	[°]	97	97	97
Y	Seat tilt angle	[°]	15, 20, 25	15, 20, 25	15, 20, 25
Z	Footrests tilt angle	[°]	95-150	95-150	95-135
	Maximum user weight	[kg]	40	60	90

### WEIGHT AND DIMENSIONS

Measurement		Size 1	Size 2	Size 3
Width	[cm]	63	68	73
Width after folding	[cm]	63	68	73
Length	[cm]	112	126	152
Length after folding	[cm]	87	91	100
Height	[cm]	97	109	123
Height after folding	[cm]	48	48	54
Weight	[kg]	18,8	21	26,6

**Note:** Measurement tolerance range is +/- 1 cm / ½ in..

**Note:** Approximate values, depending on the product configuration.

\*The value cannot be main factor for selection of product size





## 5. HOW TO USE THE STROLLER

### KEY



**Info icon** – appears in the manual when you need to read the additional information.



**Info icon** – appears in the manual when you need to press the button shown in the photo.



**Info icon** – appears in the manual when you need to use the tools indicated in the assembly description.



**Warning icon** – appears in the manual when you need to take special precautions.





## 5.1. FOLDING THE STROLLER



Keep the user away when folding and unfolding the stroller as there is a serious risk of injury.



Remove the safety rail before folding the stroller.






You can disassemble the wheels, e.g. for transportation (see chapter 5.4.1, page 16; 5.4.2, page 17).

## 5.2. UNFOLDING THE STROLLER





Make sure that the stroller has been correctly unfolded, the elements of the folding mechanism overlap the frame and are blocked. You should hear the characteristic "click" sound when the elements are overlapping.



When unfolding the stroller, make sure that the upholstery does not get between the elements of the folding mechanism!



## 5.3. FOOTREST

### 5.3.1 FOOTREST TILT ANGLE ADJUSTMENT



At the lowest footrest position, a collision with the front wheels may occur.

### 5.3.2 FOOTREST LENGTH ADJUSTMENT



#### REQUIRED TOOLS

1 hex key 3





### 5.3.3 FOOTPLATE TILT ANGLE ADJUSTMENT



#### REQUIRED TOOLS

- 1 spanner key 13 



### 5.4. WHEELS

#### 5.4.1 REMOVING REAR INFLATABLE WHEELS



In order to remove the rear inflatable wheels, push the button located in the center of the wheel.





## 5.4.2 REMOVING FRONT SWIVEL PU WHEELS



### REQUIRED TOOLS

**1** hex key 5





5.4.3 MOUNTING FRONT INFLATABLE WHEELS / SMARTRACTION SYSTEM



REQUIRED TOOLS

1 spanner key 17 

2 spanner key 22 



25



26




27




28



29



30



31



**SmarTraction System** – It is a system which makes it possible to change the position of the wheel axles depending on the user's needs. You can move them forward or backward due to special holes in the stroller frame, where you put the bushing, into which the wheel axles are inserted and fastened in the frame with the use of the quick-release system. Thanks to this solution, it is possible to adjust the distance between the front and rear wheel axles in an easy and simple way. There are two modes to choose from: **CITY MODE** and **CROSS MODE**. **The change of the wheel axles position is possible only in case of the rear inflatable wheels. It is not possible in case of the front swivel PU wheels.**



**CITY MODE**



**The stroller is more maneuverable and mobile**  
By mounting the front and rear wheel axles as close as possible, you can achieve a greater maneuverability and mobility of the stroller. This configuration of wheel axles is recommended on an even and flat asphalt road.



**CROSS MODE**

**The stroller is more stable**

By mounting the front and rear wheel axles as far as possible, you can achieve a greater stability of the stroller and a smaller risk that the stroller will fall over. This configuration of wheel axles is recommended on an uneven and bumpy terrain.





When changing the mounting place of the wheels in the stroller, the distance between its axles is also changed. Remember that after changing the mounting place of the wheels, they should be arranged symmetrically on the right and left side of the stroller.

#### 5.4. CENTRAL BRAKE



**34**



**35 WHEELS BLOCKED**



**36 WHEELS UNBLOCKED**



The central brake only works when the inflatable wheels are mounted on the stroller. It does not work with the indoor TANGO swivel wheels mounted in the HOME MODULE.



**37**

## 5.5. ADJUSTABLE BACKREST



Do not adjust the backrest with the user in the stroller.



There are three backrest positions possible.





There is a pocket at the rear of the backrest, where you can keep small, necessary accessories. Do not overload the pocket as it may cause the stroller to lose its balance. The maximum load of the pocket is 1,5 kg.



43

## 5.6. TILT-IN-SPACE FUNCTION



44




45



46



Due to the fact that the seat of the stroller is connected to the backrest, while adjusting the backrest, the tilt angle of the seat is also changed (tilt-in-space function).




## 5.7. HEADREST






## 5.8. 5-POINT SAFETY BELTS







The 5-point safety belts can be mounted at the desired height, depending on the user's needs.







 The strap should be threaded through the sliding buckle as shown in the photo. The length of the strap depends on the position of the sliding buckle.



 Always make sure that the user has the 5-point belts properly fastened while using the stroller. In order to check that, pull the belts firmly in different directions and see if the belts do not unfasten.

## 5.9. UPHOLSTERY

### 5.9.1. FLEX UPHOLSTERY



The flex upholstery is the supporting element of the stroller. Each time before use, control its condition and check if it is properly fastened and does not unfasten under pressure. If the flex upholstery has any sign of damage or loosening of stitching, contact the manufacturer to have the product repaired. Only an undamaged and properly fastened upholstery guarantees safety when using the stroller.











The flex upholstery is the supporting element of the stroller and has 3 points which hold up the structure (fig. 95). Each time you mount the flex upholstery on the stroller, check if it is properly fastened





### 5.9.2. WATERPROOF UPHOLSTERY HYDRO™





Both the flex and waterproof HYDRO™ upholstery can be cleaned with water and mild cleaning products. **Remember to remove the wooden element from the inside of the upholstery before cleaning!**

## 5.10. SAFETY RAIL

i

The front rail has a safety lock which protects it against an accidental removal.



i

One end of the safety rail is shorter than the other, which makes it possible to open the rail easily and quickly.





The height of the safety rail is adjustable. To change the height of the safety rail, place the screw in the hole at your choice. There are 5 height levels available.



### REQUIRED TOOLS

1 hex key 3



## 6. ACCESSORIES

### 6.1 NVE\_627 NOVA MODULE



In order to convert the NOVA EVO stroller into the NOVA stroller, mount two inflatable wheels at the front of the stroller (1 x NOVA MODULE). The description how to mount the inflatable wheels is available in the chapter 5.4.3 of the user manual (**see chapter 5.4.3, page 18**).

Before mounting the inflatable wheels, remove the front PU swivel wheels (**see chapter 5.4.2, page 18**).

## 6.2 NVE\_625 HOME MODULE

### HOME MODULE



In order to convert the stroller into the NOVA HOME stroller, mount 4 indoor TANGO swivel wheels (**1 x HOME MODULE**).

Before mounting the indoor TANGO swivel wheels remove the rear inflatable wheels and the front PU swivel wheels (**see chapter 5.4.2, page 17**).

**MOUNTING INDOOR TANGO SWIVEL WHEELS (HOME MODULE)**

**REQUIRED TOOLS**

**1** hex key 5 



Make sure to mount the indoor TANGO swivel wheels with a brake at the front, and indoor TANGO swivel wheels with the fixed direction of movement and a brake at the rear of the stroller.



Front indoor TANGO swivel wheels with a brake.



Rear indoor TANGO swivel wheels with the fixed direction of movement and a brake.






At the lowest footrest position, a collision with the front wheels may occur.

### 6.3 NVE\_128 ABDUCTION BLOCK



### 6.4 NVE\_405 FOLDING CANOPY





### 6.5. NVE\_445 FOOTMUFF




**136**

**137**

## 6.6. NVE\_114 H HARNESS NVE\_130 4 POINTS SAFETY VEST

**i**

Before mounting the H Harness or the 4-point safety vest remove the standard 5-point safety belts (see chapter 5.8, page 24, in reverse order) as well as the waterproof HYDRO™ upholstery (see chapter 5.9.2, page 31, in reverse order).


**138**

**139**


**140**

**141**




  
**Make sure that the belt is placed behind the wooden element.**


  
The waterproof HYDRO™ upholstery should be put back on the stroller as described in the user manual (see chapter 5.9.2, pages 31-32).





Unbuckle the bottom strap of the vest. Remove both the strap buckle and the sliding buckle (fig. 150-152). After removing the buckles, mount the strap as shown in the photos below fig. 153-160).




**154**


**155**


Put the sliding buckle (fig. 156) and the strap buckle (fig. 157-158) on the strap again. Insert the end of the strap into the sliding buckle (fig. 159-160). Fasten the strap buckle with the bottom part of the vest (fig. 161-162).


**156**

**157**

**158**

**159**

**160**



Mount the upper straps of the vest on the stroller in the same way as the 5-point safety belts are mounted (see chapter 5.8, pages 24-25, fig. 55-63).



**Before using the vest, read the safety rules included in the user manual. Information about the safe use of the vest can be found on page 47 of the user manual. Improper use of the vest and ignoring the safety rules may endanger the user's health and life.**



## 7. THE USE OF STROLLER



**NOVA EVO®** stroller has to be used with compliance to guidelines included in this user manual. Remember to use the product only with someone who is trained and knowledgeable. Do not forget to prepare both product as well as patient before using stroller!

### 7.1. PREPARATION OF THE PATIENT

**Before using the product, it is essential to make sure that:**

- patient's outfit is light and comfortable sporty clothes are advisable;
- shoes should be comfortable and, moreover, they should stabilize ankle;
- due to hygienic reason and to avoid abrasion it is advisable for the patient to always wear clothes when using the device, skin must not be naked especially in places of regular contact with upholstery;
- drains, bags, catheters etc. if applied, have to be secured to function properly and not to disturb during using the product;
- if it is necessary for a patient to wear orthopedic supply (collar, corset, orthosis, shoes) if it is recommended by doctor.



Outfit should protect patient's skin from direct contact with upholstery as well as other parts of equipment.

Only doctor can make a decision concerning using any kind of orthopedic supply!



## 7.2. PREPARATION OF THE STROLLER

### 1. Before use it is necessary to:

- Unfold the device;
  - Place the device on even ground;
  - Make sure that all parts are serviceable, correctly mounted and stable;
  - Block brakes;
  - Adjust the device to patient's measurements;
  - Undo safety belts or vest so as to place easily the patient in the device;
  - Place patient in the device;
  - Fasten all elements which stabilize patient e.g. vest, safety belts
  - Mount additional equipment e.g. abduction block, tray.
2. In case of patients who have problem with bladder, it is necessary to put special cloth on seat. This cloth can be purchased in pharmacy or medical shops.



In case of any obstacles such as step or curb, it is necessary to raise front wheels by pressing the handle down. Next, place front wheels on step and then raise back part of stroller with the aid of handle. Be careful not to tilt the stroller back too much. Remember to place stroller forwards.

It is essential to hold stroller with the aid of handle during using it otherwise stroller can roll.

Before the product is used, it is highly necessary to make sure that all belts as well as safety equipment are correctly fastened as well as adjusted. What is more, check if they do not become undone due to the pressure and, moreover, make sure that clampings from adjustable elements are firmly tightened!

Equipment works correctly only when the device is adjusted according to patient's height as well as weight. Remember to take all these measurements into account while purchasing the product. Correct adjustment is the most important due to the fact that only then the product is safe. It is forbidden to force child to keep specific position because he/she will feel uncomfortable and it may lead to wrong posture.



## 7.3. HOW TO USE VEST



The vest provides posture support for wheelchair users as well as patients sitting in positioning chairs or other rehabilitation seats. When properly used, the vest allows the patient to adopt the correct sitting position for a longer period of time. The design and shape of the vest make it possible for the patient to be free to move their upper limbs within the constraints of the patient. **If the vest is not used properly, it can pose a serious threat to the health and life of the patient!**

### Possible threats incident to the use of the vest:

- 1) **If the lower vest straps are fastened too loosely, there is a risk that the vest will move upwards to the patient's neck, making it difficult to breathe** – always check the tension of the straps and set them in such a way that the lower edge of the vest is located near the lower abdomen of the patient.
- 2) **In a situation when the patient slips off the seat, the upper edge of the vest covers the patient's neck and blocks the trachea, making it difficult to breathe** – it is therefore necessary to ALWAYS use the accessories stabilizing the patient's pelvis, such as a pelvic belt, thigh abduction belts or an abduction block.
- 3) **There is a risk that the upper straps of the vest will press the patient's neck, blocking the trachea and making it difficult to breathe.** This might happen if the patient's body is improperly secured in the seat or stroller – it is then necessary to use the optional accessories of these devices and stabilize the patient's body using the lateral supports. If there is no possibility to use such accessories with the device and the patient's body cannot be properly stabilized, the device should be replaced with the one having such an option.
- 4) **The vest is fastened too tightly to the backrest and as a result the patient cannot breathe easily** – in this case the vest should be loosened in such a way that it is possible for the caregiver to insert a straightened hand between the vest and the patient's body. If, however, we want to stabilize the patient's body properly, we should use a dynamic orthopedic vest, e.g. a Flex-point vest, which does not cause the breathing difficulty for the patient.

**Instructions for the proper use of the vest:**

- 1) Make sure that the patient is sitting in the best possible position:
  - in the middle of the seat,
  - deeply in the seat in such a way that the patient's lumbar spine rests against the backrest.
- 2) Secure the patient's pelvis with, for example, a pelvic belt, thigh abduction belts or an abduction block.
- 3) Make sure that the side supports of the device, if used, are set in the correct position. It is required that the patients with the poor trunk control use the side supports. The side supports limit the lateral movements of the patient's body, which could lead to the compression of the vest straps on the neck of the patient.
- 4) Secure the vest:
  - fasten and adjust the lower straps of the vest, so that the lower edge of the vest is located near the lower abdomen of the patient,
  - fasten the horizontal strap of the vest (if present),
  - fasten the shoulder straps of the vest,
  - check if the buckles of the vest do not open under pressure.
- 5) Check the patient's position:
  - make sure that the vest cannot move by itself towards the patient's neck,
  - make sure that you can insert a straightened hand between the vest and the patient's body.

**To prevent injuries and to minimize the risks associated with the use of the vest, the following rules should be followed:**

- 1) The position of the vest **MUST** be checked each time after putting the patient in the seat.
- 2) Special attention should be paid to protecting the patient from falling forward, slipping off the seat or lifting the pelvis. The vest **MUST** be used with the accessories which will secure the patient's pelvis at all times (e.g. the pelvic belt, thigh abduction belts or the abduction block). The use of the vest cannot be considered separately but as a part of a general posture support system including supports such as pelvic belts, foot support, torso stabilization, thigh abduction belts, tilt-in-space position and reclining backrest.
- 3) Special attention should be paid to ensuring that the entire foot of the patient has stable support. If the device is not equipped with a foot support, the patient's feet must rest firmly against the ground.



- 4) It should be remembered that any change of the inclination angle of the backrest, seat, foot support or tilt-in-space function might affect the effectiveness of the belts and supports!
- 5) Special attention should be paid to the use of stabilization elements or side supports. If the patient has poor trunk control, the side supports will limit the lateral movements, which could lead to the compression of the vest straps on the neck of the patient.
- 6) The mounting of the vest and adjusting it to the patient can only be carried out by a suitably trained person, who is familiar with the instructions for the use of the vest and the rules of its safe use.
- 7) Do not make any repairs or modifications to the vest on your own. In order to repair or replace any part of the vest, please contact the manufacturer.
- 8) Do not use a damaged or incomplete vest, as this poses a serious threat to the health and life of the patient!
- 9) Do not leave the patient in the device unattended!
- 10) The correctness of the use of the vest should be particularly controlled in accordance with the recommendations in case of children, patients with impaired memory or mentally ill.
- 11) Make sure that all people taking care of the patient are familiar with the rules of the safe use of the vest.

### **How to care for the vest:**

The vest works effectively only if kept clean and maintained in accordance with the manufacturer's instructions. Before each use of the vest, check whether:

- 1) The vest straps are not frayed or worn. If so, have them exchanged for new ones,
- 2) The vest straps are not attached to the wheels of the device,
- 3) The buckles are not damaged. If so, have them exchanged for new ones,
- 4) The Velcro fastening holds firmly. If not, remove any dirt or hair from the Velcro, or have it replaced with a new one,
- 5) There have been no changes in the patient's height or weight from the moment the vest was adjusted to the patient,
- 6) The shoulder straps of the vest are adjusted so that they are not too loose or too tight,
- 7) The vest should be cleaned regularly with the use of warm soapy water and dried out of the direct sunlight and other sources of heat. The vest must always be dry before usage.

**If you have any doubts or questions about the use of the vest, please contact your doctor or physiotherapist.**



## 8. SAFETY RULES



**NOVA EVO®** stroller has to be used in compliance with guidelines in this user manual.

It is highly important to consult either doctor or physiotherapist in order to make right decision concerning the choice of stroller and its size. It is also essential to familiarize with the guidelines described in user manual before using the device.

### 1. Before use it is necessary to:

- make sure that there are no damages or different faults which could endanger patient's life; check if all parts are working and are correctly mounted and, moreover, check if they are not broken;
- check if all screws, knobs and all glued, sewed or tapped parts are tightened;
- check if brakes are operating correctly;
- check if buckles are operating correctly;
- always use safety belts, Velcro tapes should overlap min. 8cm;
- adjust the device according to patient's measurements and needs e.g. the adjustment of pelottes, footplate etc.;
- make sure that all movable parts and knobs for adjustments are stable and locked
- secure patient with the aid of elements such as vest, foot stabilizing belts, pelottes, abduction block, headrest. It is highly important to secure patient because it safeguards him/her from falling out. If abduction contracture occurs, it is recommended to use abduction block in order to prevent ulcers between knees.





2. It is forbidden to use stroller for different purposes than it is intended for.
3. It is forbidden to use stroller when it is damaged or any part is missing.
4. It is forbidden to leave the patient unattended.
5. Stroller should be used only when there is someone who is trained and knowledgeable.
6. It is forbidden to leave stroller and the patient on unstable ground unattended, even if the brake is engaged.
7. It is forbidden to either go up or down the stairs; it doesn't matter if the patient is placed in the device or not.
8. It is important to engage brake during getting into and getting out of the stroller.
9. Make sure that the user does not stand on the footrest when getting in and out of the stroller.
10. Make sure that all limbs are not close to spokes or other adjustable elements.
11. The stroller should be stored and used away from sources of fire or high temperature. Moreover, stroller shouldn't be stored in rooms with high humidity.
12. Secure stroller during transportation in order to avoid any possible damages. Manufacturer doesn't take responsibility for damages which appeared due to improper preparation of the device during transportation.
13. Before using the stroller, remove all paper and plastic packaging. The packaging is not a toy! Keep away from children and animals





**It is forbidden to place the patient in the device while being in car, bus or train.**

**The product contains small parts which children may swallow.**

**Do not use a product that has any signs of damage, tearing or wear. Only an undamaged product and proper use guarantee its right and safe working. If the product is damaged, discontinue use and contact the manufacturer.**

**Under no circumstances should an unqualified person repair or alter the product and its construction, as the health and safety of the user depend on it as well as it will lead to the loss of warranty.**

**Remember that all the loads hung on the push handle may cause the stroller to lose its balance.**

**The weight of user should not be bigger than maximal weight described on page 10 where you will find information regarding this issue.**

## 9. TERMS OF WARRANTY

- 1) Akces-MED Sp. z. o. o sells products free from manufacturing defects or faults.
- 2) Manufacturer gives the purchaser **24-month** guarantee from the date of purchase. All defects revealed during warranty period will be repaired free of charge.
- 3) All defects which appear after warranty period can be repaired only when the user will pay for it. This rule doesn't apply to refunded products - user can ask National Health Refund for refund.
- 3a) The stroller is covered by a 36-month warranty on the territory of Canada – the details of the rights resulting from the extended warranty are regulated by the warranty card.
- 4) Manufacturer claims that the product can be used for **48-months**. User can use this product after this period only when he/she will send the device to manufacturer so as to service it. Next, manufacturer determines again how long the product can be used and, moreover, the date of the next service will be determined.
- 5) If product is damaged, please contact our dealer or Akces-MED Sp. z. o. o customer service (ph. 17/864 04 70, e-mail: [export@akces-med.com](mailto:export@akces-med.com))





## 10. MAINTENANCE AND CLEANING

The stroller **NOVA EVO**®, which is a medical device, should always be kept clean and should be used in accordance with the manufacturer's instructions. The stroller should be kept in a clean, dry and ventilated room. Maintenance and cleaning must be done at least once a month. Periodic cleaning and maintenance of the chair is essential to keep its proper functioning and safety.

### **Recommendations for cleaning the upholstery and other elements of the stroller:**

1. Wash with warm water using a soft cloth and mild cleaners.
2. Make sure the cloth is not too wet, as this may enlarge the stain and damage the upholstery or foam inside the upholstery.
3. Remember to remove excess moisture from the upholstery before the next cleaning.
4. After cleaning the stains, the upholstery should be rinsed with lukewarm water and dried with a clean cloth.
5. You must not soak the upholstery.
6. The upholstery must not get wet during cleaning.
7. Do not use hard brushes for cleaning the upholstery.
8. Do not clean the upholstery with strong chemicals, e.g. solvents or chlorine-based caustics.
9. Do not clean, replace any elements or maintain the chair during use.
10. Metal and plastic elements should be cleaned by hand with suitable cleaners, after having read the recommendations of the manufacturer of the given cleaner.



**The manufacturer of the stroller, Akces-MED Sp. z o. o., is not responsible for any damages caused by using unsuitable cleaners.**





Remember that upholstery has to be always dry before using!

**Before using it is necessary to:**

- clean wheels;
- check if base as well as wheels are correctly mounted; check if brake stops turning wheel;
- check if wheelchair construction is free of breaks or any kinds of deformities;
- check if all movable parts are working and all screws, nuts, knobs etc. are correctly tightened;
- check if belts, vests, trunk belts, feet stabilizing belts etc are not frayed; check if buckles and clamps work correctly.



Manufacturer takes responsibility only for hidden faults or damages which appeared despite the fact that user complied with the guidelines.

In case of any damages, it is essential to contact the person who sold you the product.





## 11. LABELLING

### 11.1. SYMBOLS

 Declaration that the product meets the essential requirements

 Catalog number

 Serial number

 Please read the manual instructions before use

 Do not iron

 The product does not contain latex.

 Store in a dry place

 Careful! delicate content

 Protect from sunlight

 Storage temperature limitation

 Relative humidity storage

 Do not chlorinate

 Do not dry in a drum dryer

 Do not bleach

 Dry vertically

 Prohibition of use without fastened safety belts

 Prohibition to leave the patient without care during use of the product

 It is forbidden to either go up or down the stairs; it doesn't matter if the patient is placed in the device or not.

 Manufacturer

 Production date

 You need to read the additional information.

 You need to press the button shown in the photo.

 You need to use the tools indicated in the assembly description.


 You need to take special precautions.

## MANUFACTURER



**Akces-MED Sp. z o.o.**

Jasionka 955B  
36-002 Jasionka  
tel. (17) 864 04 70  
e-mail: [export@akces-med.com](mailto:export@akces-med.com)  
BDO: 000031619



Authorised Dealer



**akces-med.com**

Date of issue: 09.12.2019  
The last update: 30.03.2020  
© Sławomir Wroński 2020