



**HIPPO+™**  
special stroller



Special needs stroller **HIPPO+™** is equipped with many interesting and useful solutions. It has been adapted to changing weather conditions thus it will protect the user from sunlight, wind, or snow depending on the season. Perfect for long tours.

Modern design of the stroller makes it look very elegant, and yet it still retains all its therapeutic functions. Stylish appearance of **HIPPO+™** does not even resemble a special needs stroller.

One of its main advantages is a reversible seat. As a result of this solution, the patient is able to observe the surroundings (in the case of forward facing) or caregiver (in the case of rear facing).

The product comes standard with many useful accessories such as capacious basket or cover for feet. The whole is made of the finest materials that ensure both durability and reliability.

Available in 2 sizes. Especially recommended for youngest.



### IT IS RECOMMENDED FOR CHILDREN IN THE FOLLOWING CASES:

- Infantile cerebral palsy – different forms.
- Meningomyelocele.
- Brain and spinal cord injuries.
- Muscle diseases connected with paresis and paralysis.
- Genetic disorders connected with paresis.
- Other diseases connected with paresis, paralysis and injuries regarding locomotor system.
- After traumatic brain injury, spine with spinal cord and lower limbs injury.



### STANDARD EQUIPMENT:

Foldable frame with wheels and rear brake | Swivel front wheels 360° | Drive system with shock absorbers | Under seat storage basket | Abduction block | Reclining Backrest | Adjustable side supports | 5 points seat belts | Adjustable headrest | Adjustable and folding canopy | Detachable and reversible seat unit (traveling face to face with the child) | Adjustable push handle made of ECO leather | Leg cover | Pumped wheels front/rear | Front rail | Rain cover | Bag | Pump | Reduction set consisting of additional seat (size 1) | Lumbar pad (size 1)



ADVANTAGES:



**Stable structure**

The product ensures adequate stability and safety of a Patient.



**Movable seat**

You can turn the seat easily so that the child may be facing towards caregiver or movement direction.



**Lightweight construction**

The product is made of light metal alloys which positively affects the low weight of the device.



**Easy to use**

Intuitive way of adjustment and easy assembly.



**Reflective elements**

The product has reflective elements to ensure your safety.



**24 months warranty**

We offer 24 months warranty on the product. We also provide after sales service.



Front swivel wheels with directional lock facilitate easy maneuvering of the stroller.



Warm cover works great during bad weather conditions.



Stylish and capacious bag can store many useful little things.

COLORS



Mint



Anemon



Chilli



Lotus



*Our passion - your independence!*

Colors may differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory. Current offer on our website.

Akces-MED sp. z o.o.  
Jasionka 955B  
36-002 JASIONKA

export@akces-med.com  
+48 (17) 864 04 77  
www.akces-med.com



Symbol	Measurement	Unit	Size 1	Size 2
A	Seat width	[cm]	20*/ 31	34
L	Seat depth	[cm]	28	32
	Backrest height	[cm]	40*, 45	55
L1	Seat depth to the abduction block	[cm]	18	22
C	Headrest height	[cm]	40-51*/ 44-54	55 - 64
E	Backrest width	[cm]	31	34
F	Footrests length	[cm]	20*, 24*/16, 20	22, 26, 30
	Footer length	[cm]	14	18,5
X	Backrest tilt angle	[°]	90-140	90-140
Y	Seat tilt angle	[°]	10	10
Z	Footplate tilt angle	[°]	90-180	90-180
	Max. user weight	[kg]	20	30
	Indicative user height**	[cm]	70-100	90-120
Product dimensions & weight			Size 1	Size 2
	Width	[cm]	66	66
	Length	[cm]	109	113
	Height	[cm]	105	113
	Width after folding	[cm]	66	66
	Length after folding	[cm]	115	125
	Height after folding	[cm]	56	56
	Weight	[kg]	19,5	22

**Note:** Measurement tolerance range is +/- 1 cm / 1/2 in.  
 \*Measurement starts from narrowing pad.  
 \*\* The value cannot be main factor for selection of product size.




## Safety



### Led lights (2 pieces)

Front and rear LED lights for increased safety.



### Car seat

Designed to carry a child in the car. It provides a safe and comfortable journey.



## Stabilization



### Head stabilizer with belt

Mounted in place of the standard headrest. It can be precisely adjusted to the user's needs. Soft head belt additionally stabilizes the forehead.



### Cushion to narrow down the seat (size 2)

The cushion narrows down the seat. It is made of elastic foam trimmed with soft material. 3 cm wide



### Lumbar support (size 2)

It allows better stabilization of the lumbar back.



### Thigh abduction belts

They can be easily attached to the chair by use of Velcro straps. Their use allows to keep legs in abduction and also prevents the patient from slipping off of the device.



### H vest

It acts as a trunk stabilizer. HOLD&PULL© system makes it easier to adjust the vest to the patient body. You can quickly detach the vest if necessary.



### 4 point vest

It acts as a trunk stabilizer. Wide range of adjustment helps to fit it to the user's needs. Equipped with HOLD&PULL© system.



### Depth reducing pad (3 cm / 5 cm)

It allows to reduce the depth of the seat. Made of elastic foam material. Available widths of 3 cm and 5 cm.



## Comfort




### Umbrella

Attachable to the stroller's frame, perfectly protects from the sunlight or rain.



### Cup holder

An universal holder mounted to the frame. It enables easy and quick access to a bottle of drink.



### Tray made of plastic

Perfect for working and taking meals, it also supports the therapy.



### Gondola

Designed to carry the child in a lying position. It can be mounted to the frame without additional adapters.



**Headrest cover**

It helps to keep the headrest clean. Easy to assemble/disassemble, it can be washed.



**Wheel covers (4 pcs.)**

Use it on the dirt tires when you want to take the products indoor or transport in the trunk. .



**Mosquito Net**

Excellent protection against insects



**Build-up**



**Kid-sit buggypod IO**

For convenient and safe transportation of a second child without having to buy a double stroller. The buggypod IO has been designed so that within a few seconds, it can be folded, unfolded, or completely disassembled. 5 point belts seat and folding footrest provide the child with comfort and safety. Using the conversion kit, you are able to convert the kid-sit into an independent three-wheeled stroller. Load capacity up to 20kg.



**Buggypod LITE**

Starters are very helpful if you want to take two children for a walk - in this case, you just need to mount the attachment to a larger trolley. Thanks to the special buckles, it can be easily adjusted to Ursus, Ombrello, Hippo+



**Buggypod PERLE adapter**

Perle is a well-designed seat that fits into the BuggyBoard Maxi.



**Buggypod Perle seat**

Perle is a well-designed seat that fits into the BuggyBoard Maxi.



**Rain cover for the kid-sits**

An excellent protection from the rain, snow, and wind either for the user and stroller



**Extension bar**

An adapter enabling to mount the Buggypod kid-sits to the strollers of other manufacturers.



**Mounting set**

It enables easy assembly of the Buggypod kit-sits to other strollers whose frame diameter is larger than 40mm.



**Conversion kit (for The buggypod kit-sits)**

Set of elements thanks to which the Buggypod kit-sit can be converted into a fully functional stroller.



Spare parts



**Rear wheel (1 piece)**

Rear wheel 12" ½ x 2" ¼. There is a possibility to order only an inner tube or tire.



**Front wheel (1 piece)**

Front wheel 10" x 2". There is a possibility to order only an inner tube or tire.



Storage



**Stroller bag**

An excellent solution to protect the folded stroller from damages (e.g. during the transportation). It perfectly fits the stroller shape and dimensions.



**Organizer**

Handy carrying bag on the handle of the trolley. Specially designed compartments make it easy to keep order by storing the most needed items at hand. The Organizer has two useful holes for a cup of coffee or a bottle of water.