



**HIPPO+**<sup>TM</sup>  
special stroller



Special needs stroller **HIPPO+™** is equipped with many interesting and useful solutions. It has been adapted to changing weather conditions thus it will protect the user from sunlight, wind, or snow depending on the season. Perfect for long tours.

Modern design of the stroller makes it look very elegant, and yet it still retains all its therapeutic functions. Stylish appearance of **HIPPO+™** does not even resemble a special needs stroller.

One of its main advantages is a reversible seat. As a result of this solution, the patient is able to observe the surroundings (in the case of forward facing) or caregiver (in the case of rear facing).

The product comes standard with many useful accessories such as capacious basket or cover for feet. The whole is made of the finest materials that ensure both durability and reliability.

Available in 2 sizes. Especially recommended for youngest.

**MEDICAL  
DEVICE**  
CLASS 1



**WYRÓB  
MEDYCZNY**  
KLASA 1



## IT IS RECOMMENDED FOR CHILDREN IN THE FOLLOWING CASES:

- Infantile cerebral palsy – different forms.
- Meningomyelocele.
- Brain and spinal cord injuries.
- Muscle diseases connected with paresis and paralysis.
- Genetic disorders connected with paresis.
- Other diseases connected with paresis, paralysis and injuries regarding locomotor system.
- After traumatic brain injury, spine with spinal cord and lower limbs injury.

## STANDARD EQUIPMENT:

Foldable frame with wheels and rear brake | Swivel front wheels 360° | Drive system with shock absorbers | Under seat storage basket | Adjustable footrest with supporting belt | Removable abduction block | Reclining backrest | Adjustable headrest | Adjustable folding canopy | Detachable and reversible seat unit | Adjustable push handle made of ECO leather | Leg cover | Pumped wheels (front/rear) | Front rail | Rain cover | Bag | Pump | Reduction set consisting of additional seat | Lumbar pad



### ADVANTAGES:



**Stable structure**

The product ensures adequate stability and safety of a Patient.



**Movable seat**

You can turn the seat easily so that the child may be facing towards caregiver or movement direction.



**Lightweight construction**

The product is made of light metal alloys which positively affects the low weight of the device.



**Easy to use**

Intuitive way of adjustment and easy assembly.



**Reflective elements**

The product has reflective elements to ensure your safety.



**24 months warranty**

We offer 24 months warranty on the product. We also provide after sales service.



Maneuvering the stroller is really easy due to the front swivel PU wheels.



A warm leg cover is perfect for harsh weather conditions.



The seat can be put in the forward or rearward position.



A proper stabilization of the user is provided by a wide range of standard equipment.

### COLORS



LOTUS



ANEMON



VULCANO



DENIM




FLAMING



AQUA

**Attention!** The lead time for selected colors (DENIM/AQUA/VULCANO/FLAMING) may be longer.



Symbol	Measurement	Unit	Size 1	Size 2
A	Seat width	[cm]	20* / 31	34
L	Seat depth	[cm]	28	32
	Backrest height	[cm]	40*, 45	55
L1	Seat depth to the abduction block	[cm]	18	22
C	Headrest height	[cm]	40 – 51*/ 44 - 54	55 – 64
E	Backrest width	[cm]	31	34
F	Footrests length	[cm]	20*, 24*/16, 20	22, 26, 30
	Footer length	[cm]	14	18,5
X	Backrest tilt angle	[°]	90-140	90-140
Y	Seat tilt angle	[°]	10	10
Z	Footplate tilt angle	[°]	90-180	90-180
	Max. user weight	[kg]	20	30
	Indicative user height**	[cm]	70-100	90-120
Product dimensions & weight			Size 1	Size 2
	Width	[cm]	66	66
	Length	[cm]	109	113
	Height	[cm]	105	113
	Width after folding	[cm]	66	66
	Length after folding	[cm]	115	125
	Height after folding	[cm]	56	56
	Weight	[kg]	19,5	22

**Note:** Measurement tolerance range is +/- 1 cm / 1/2 in.  
 \*Measurement starts from narrowing pad.  
 \*\*Height cannot be main factor for selection of product size.



Safety



**HPO\_001 LED lights**  
Front and rear LED lights for increased safety.



Positioning



**HPO\_101 Thigh abduction belts**  
They can be easily attached to the chair by use of Velcro straps. Their use allows to keep legs in abduction and also prevents the patient from slipping off of the device.



**HPO\_114 H Harness**  
It acts as a trunk stabilizer. HOLD&PULL™ system makes it easier to adjust the vest to the patient body. You can quickly detach the vest if necessary.



**HPO\_124 Lumbar pad**  
It allows better stabilization of the lumbar back.



**HPO\_130 4 points safety vest**  
It acts as a trunk stabilizer. Wide range of adjustment helps to fit it to the user's needs. Equipped with HOLD&PULL™ system.



**HPO\_137 Reduction set consisting of additional seat**  
The cushion narrows down the seat. It is made of elastic foam trimmed with soft material.



**HPO\_138 Headrest with belt**  
Mounted in place of the standard headrest. It can be precisely adjusted to the user's needs. Soft head belt additionally stabilizes the forehead.



**HPO\_151 Depth reducing pad 3 cm**  
It allows to reduce the depth of the seat. Made of elastic foam material.



**HPO\_152 Depth reducing pad 5 cm**  
It allows to reduce the depth of the seat. Made of elastic foam material.



**HPO\_161 Sandals stabilizing feet and ankles**  
They protect and stabilize the ankle joint. They are made of elastic material, recommended for users with increased muscle tone. They can be used with or without shoes.



**HPO\_162 Dynamic ankle stabilizers (FP-22)**  
Dynamic FlexPoint® ankle stabilizers provide great support when stabilizing the user's feet. Incorrect positioning of the feet, despite appropriate stabilization of the pelvis and chest, makes it impossible to keep the correct sitting position.



Comfort



**HPO\_402 Umbrella**  
Attachable to the stroller's frame, perfectly protects from the sunlight or rain.



**HPO\_409 Mosquito Net**  
Excellent protection against insects.



**HPO\_411 Cup holder**  
An universal holder mounted to the frame. It enables easy and quick access to a bottle of drink.



**HPO\_414 Tray Plexiglas**  
Perfect for working and taking meals, it also supports the therapy.



**HPO\_415 Wheel cover**  
Use it on the dirt tires when you want to take the products indoor or transport in the trunk.



**HPO\_429 Soft footrest cover**  
Comfortable pad, protecting the feet from direct contact with the footplate



**HPO\_445 Front wheel cover**  
Use it on the dirt tires when you want to take the products indoor or transport in the trunk.



**KDO/KDH\_448 Winter muffs**  
Comfortable and warm muffs protect your hands against a chill during walks.



**KDO/KDH\_448 Lap tray**  
A lap tray is a multifunctional and portable tray. A special bag filled with granules makes it possible to stably place it on your knees or on other uneven surfaces.



Storage



**HPO\_506 Travel Bag for stroller**  
An excellent solution to protect the folded stroller from damages (e.g. during the transportation). It perfectly fits the stroller shape and dimensions.



**HPO\_509 Organizer bag**  
Handy carrying bag on the handle of the trolley. Specially designed compartments make it easy to keep order by storing the most needed items at hand. The Organizer has two useful holes for a cup of coffee or a bottle of water.