



NONO™
bath seat



NONO™

bath deckchair

AKCES MED

NONO™ bath seat for children was designed to care for the child while providing the appropriate level of stability.

Thanks to its simple design, the product is very helpful when bathing, easily fits into a standard bathtub and allows for optimal positioning of a child.

Secured and comfortable position makes bathing easier.

Lightweight, aluminium frame and net upholstery protect against adverse effects of water.

All the products for bathing have pipes filled with polyurethane foam, which prevents water ingress. After removing the product from bath or shower water will not leak out which will remain your floor dry.

Seat is equipped with backrest tilt's angle adjustment which makes it possible to adjust it to child's needs (for both sitting and lying position).

Available in **2** sizes: **MINI** and **MAXI**



IT IS RECOMMENDED IN THE FOLLOWING CASES:

- Infantile cerebral palsy,
- Paralysis,
- Injury or amputation of limbs,
- Deformation of limbs,
- Degeneration or injury of joints.

COLORS



Nono MINI



Nono MAXI

STANDARD EQUIPMENT:

Waterproof upholstery | Water-resistant adjustable chair's frame | equipped with foam covers to avoid scratches on the oath tube | Pelvic belt

our passion - your independence!

Colors may be differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory.
Current offer on our website.

Akces-MED sp. z o.o.
Jasionka 955B
36-002 JASIONKA

export@akces-med.com
+48 (17) 864 04 77
www.akces-med.com

Version: 09/2017



NONO™

bath deckchair

AKCES MED®

ADVANTAGES:



Waterproof

This construction is neutral for harmful water influence.



Easy to use

Adjusting does not require the use of tools.

Lightweight construction

The product is made of lightweight and durable materials

Compact

When folded occupies very little space.

Safety

Seat belts protect the child while bathing.


24

24 months warranty

We offer 24 months warranty on the product. We also provide after sales service.



Innovative articulated system with a wide range of frame adjustments makes it possible to adjust the seat to patient's measurements as well as needs.



Water-resistant frame.



Trunk support with Velcro straps.



our passion - your independence!

Colors may differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory.
Current offer on our website.

Akces-MED sp. z o.o.
Jasionka 955B
36-002 JASIONKA

export@akces-med.com
+48 (17) 864 04 77
www.akces-med.com

Version: 09/2017



NONO™

size table

AKCES MED®



Symbol	Measurement	Unit	Mini	Maxi
A	Product width	[cm]	37,5	37,5
B	Seat depth	[cm]	30	30
C1	Backrest height	[cm]	45	58
	Max user weight	[kg]	30	30

Product dimensions and weight		Mini	Maxi
Width	[cm]	40	40
Width after folding	[cm]	40	40
Length	[cm]	65	65
Length after folding	[cm]	73	84
Height	[cm]	49	59
Height after folding	[cm]	14	14
Weight	[kg]	2	2,5

Note: Measurement tolerance range is +/- 1 cm / ½ in.

our passion - your independence!

Colors may be differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory.
Current offer on our website.

Akces-MED sp. z o.o.
Jasionka 955B
36-002 JASIONKA

export@akces-med.com
+48 (17) 864 04 77
www.akces-med.com

Version: 09/2017