



NONO™
bath seat



NONO™ bath seat for children was designed to care for the child while providing the appropriate level of stability.

Thanks to its simple design, the product is very helpful when bathing, easily fits into a standard bathtub and allows for optimal positioning of a child.


Secured and comfortable position makes bathing easier.

Lightweight, aluminium frame and net upholstery protect against adverse effects of water.

All the products for bathing have pipes filled with polyurethane foam, which prevents water ingress. After removing the product from bath or shower water will not leak out which will remain your floor dry.

Seat is equipped with backrest tilt's angle adjustment which makes it possible to adjust it to child's needs (for both sitting and lying position).

Available in **2** sizes: **MINI** and **MAXI**



IT IS RECOMMENDED IN THE FOLLOWING CASES:

- Infantile cerebral palsy,
- Paralysis,
- Injury or amputation of limbs,
- Deformation of limbs,
- Degeneration or injury of joints.

COLORS



Nono MINI



Nono MAXI

STANDARD EQUIPMENT:

WATERPROOF UPHOLSTERY HYDRO™ | Water-resistant adjustable chair's frame | equipped with foam covers to avoid scratches on the bath tube | Pelvic belt



ADVANTAGES:



Waterproof

This construction is neutral for harmful water influence.



Lightweight construction

The product is made of lightweight and durable materials



Compact

When folded occupies very little space.



Easy to use

Adjusting does not require the use of tools.




Safety

Seat belts protect the child while bathing.




24 months warranty

We offer 24 months warranty on the product. We also provide after sales service.



Innovative articulated system with a wide range of frame adjustments makes it possible to adjust the seat to patient's measurements as well as needs.




Water-resistant frame.



Trunk support with Velcro straps.





Symbol	Measurement	Unit	Mini	Maxi
A	Seat width	[cm]	37,5	37,5
B	Seat depth	[cm]	30	30
C ₁	Backrest height	[cm]	45	58
	Max. user weight	[kg]	60	75
	Indicative user height*	[cm]	90	110
	Product dimensions & weight		Mini	Maxi
	Width	[cm]	40	40
	Width after folding	[cm]	40	40
	Length	[cm]	65	65
	Length after folding	[cm]	73	84
	Height	[cm]	49	59
	Height after folding	[cm]	14	14
	Weight	[kg]	2	2,5

Note: Measurement tolerance range is +/- 1 cm / ½ in.