



**NOOK™**  
positioning chair



It is a corner positioning chair which helps children and young people to adopt the correct sitting position with upright legs. It is worth saying that this position enforces the correct alignment of hips thanks to the Efficient Tilt-in-space adjustment

The use in the product of an innovative four-bar linkage, supported with the gas spring, greatly improves smooth Tilt-in-space adjustment..

It is recommended as the first therapeutic chair.

Ergonomic side supports and Efficient Tilt-in-space adjustment

The use in the product of an innovative four-bar linkage, supported with the gas spring, greatly improves smooth Tilt-in-space adjustment. with pelvic belt not only make the patient feel secure but also control the correct posture. Additionally, a vast array of adjustments is beneficial due to the fact that children grow fast.

There are 3 sizes and therefore it is perfect for children as well as



### IT IS RECOMMENDED IN THE FOLLOWING CASES:

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



### STANDARD EQUIPMENT:

Base | Seat | Removable abductor | Pelvic support adjustable | Pelvic belt | Chest support adjustable | Chest belt | Castors with brakes (size 3 only)

*Our passion – your independence!*

Colors may differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory. Current offer on our website.

Akces-MED sp. z o.o.  
Jasionka 955B  
36-002 JASIONKA

export@akces-med.com  
+48 (17) 864 04 77  
www.akces-med.com

Version: 07/2016



## ADVANTAGES:



### Stable structure

The device ensures adequate stability and safety of Patient.



### Easy to use

Intuitive way of adjustment and easy assembly.



### Eco-friendly materials

We use only human and environment friendly materials.



### Spacer upholstery

The device is equipped with breathable, hypoallergenic 3D spacer fabric. Now you can be sure that the skin stays dry for longer!



### 24 months warranty

We offer 24 months warranty on the product. We also provide after sales service.



Ergonomically shaped side supports control the correct position in the chair.




The tray is perfect for different activities e.g. playing or taking meal.




Large, profiled Efficient Tilt-in-space adjustment  
The use in the product of an innovative four-bar linkage, supported






## BREATHABLE UPHOLSTERY


- Lightweight
- Hypoallergenic
- Natural protection against mold
- The antibacterial properties
- Easy to clean
- Easy decontamination
- Wide range of hardness and elasticity
- Ecological
- High durability
- Resistance to crushing




Grass  
Size 1



Sea  
Size 2




Latte  
Size 3



Red  
Optional

Optional colors available for special request. Production time ca. 14 days.



## WATERPROOF UPHOLSTERY

- Lightweight
- Easy to clean
- Soft
- Waterproof
- Easy to clean
- Ecological
- High durability
- Easy decontamination




Green Grass  
Optional



Blue Sea  
Optional




Brown Sugar  
Optional



Hot Red  
Optional

Optional colors available for special request. Production time ca. 14 days.




## MEDICAL FABRIC

- Medical Fabric
- Breathable
- Waterproof
- Fire retardant\*

Upholstery is made of innovative fabric, which is very safe for patient and easy to disinfection and because of that the product can be used by a many patients. Additional asset of the fabric is breathable and water-resistant. This is the best solution as an waterproof cover of anti bedsore mattresses in hospitals.

\*For special request



MED Blue  
Optional



MED Latte  
Optional


Optional colors available for special request. Production time ca. 14 days.

*Our passion - your independence!*

Colors may differ from the presented images. The manufacturer reserves the right to make changes in products without updating the directory. Current offer on our website.

Akces-MED sp. z o.o.  
Jasionka 955B  
36-002 JASIONKA

export@akces-med.com  
+48 (17) 864 04 77  
www.akces-med.com



Symbol	Measurement	Unit	Size 1	Size 2	Size 3
A	Seat width	[cm]	50	56	67
B	Seat depth	[cm]	29	32	41
C	Backrest height	[cm]	33-41	47-58	51-66
E	Pelvic support width	[cm]	26	30	36
E1	Chest support width	[cm]	26	30	40
K	Seat depth to the abduction block	[cm]	12-22	14-24	25-35
	Max user weight	[kg]	25	35	50
Product dimensions and weight		Unit	Size 1	Size 2	Size 3
	Width	[cm]	50	56	67
	Lenght	[cm]	60	62	72
	Height*	[cm]	38	52	58*
	Weight	[kg]	6	8	11,5

\*In case of mounted wheels dimension increases by 8 cm.

**Note:** Measurement tolerance range is +/- 1 cm / 1/2 in.



## Safety



### Soft cover for tray

Soft cover protects the user from hitting the tray. Indicated for patients with impaired limbs coordination.



## Stabilization



### Adjustable headrest

Headrest adjustable in three dimensions. Due to its special design, it can be precisely adjusted to the user's needs.



### Extended headrest

For patients with increased head circumference. It provides better side stabilization.



### Head supporting belt

Attachable to the headrest by use of Velcro straps. Used for people with limited control of the head to prevent its falling.



## Comfort



### Cotton cover for headrest

It helps to keep the headrest clean. Easy to assemble/disassemble, it can be washed.



### Swivel castors (5pieces)

Mounted to the base. They allow easy movement between rooms.



### Tray

Made of wood. It supports the therapy and allows you to play, eat, and work.



### Tray for non verbal communication

It is a very important element in the rehabilitation of patients with speech dysfunction and, moreover, it makes it easier to perform daily activities. The set consists of 38 cards that may be put in random order into 24 places of the tray. It is covered with transparent and aesthetic plate made by Plexiglas®. Parents and child can create their own symbols that they believe are significant for their communication.