



# SYSTEM URSUS™

multifunctional device



URSUS™  
special stroller



URSUS HOME™  
special stroller




URSUS ACTIVE™  
active special stroller



Ursus™ seat is a modern, safe and functional solution for patients, which prefers a high comfort and stabilization. **TILT-IQ™** function, Tilt-in-space **ClickClack™** are some of many practical solutions of **URSUS SYSTEM™**

- **Outdoor frame Ursus™** Means tours and travels in very comfortable conditions either for a child and caregivers. It provides a convenient ride in varied terrain. Swivel front wheels ensure smooth changing of direction.
- **Hi-Low frame Ursus Home™** Is a comfortable and functional solution in home conditions. Adjustable seat height greatly facilitates contact with peers, and push handles make maneuvering, even in tight spaces, much easier.



one seat  
many possibilities

- **Active frame Ursus Active™** is a mode that gives your child a freedom of movement. Thoughtful solution enables independent movement of the user or with assistance of a caregiver. Now, therapy, learning, and fun will be even more enjoyable and effective.
- **Pneumatic Hi-Low frame Quark™** Will convert **URSUS™** seat into a comfortable rehabilitation chair. Innovative system of height adjustment en-sures smooth and easy operation. Included push handle allows easy relocation of the device.

Choose  
your  
upholstery

## BREATHABLE UPHOLSTERY ATMO™



BREATHABLE



EASY TO CLEAN



DURABLE



ECOLOGICAL



ANTIALLERGIC



ANTIMICROBIAL  
& ANTIFUNGAL

## WATERPROOF UPHOLSTERY HYDRO™



WATERPROOF



EASY TO CLEAN



DURABLE



ECOLOGICAL



FIRE RETARDANT



Ultralight stroller **URSUS™** is designed for daily rehabilitation of children with locomotor disabilities. As far as its construction is concerned, the frame and seat are made of high-strength aluminum, and the chassis is additionally cushioned. **Ursus™** is equipped with **ClickClack™** feature, which means that you can easily swap the seats and adjust the stroller to the most suitable mode. The device has a wide range of adjustments such as width, seat depth, backrest height, backrest tilt angle, infinitely adjustable tilt-in-space function – **TILT-IQ™**, reversible seat. **Ursus™** has also been equipped with spatial, breathable upholstery that prevents excessive sweating as well as bedsores, and more to the point, positively influences the comfort of use. Moreover, the stroller can easily be converted into a lounger so that the child can sleep in it.

The product is available in **3** different sizes with a very wide range of accessories that will help you to meet the individual needs of user.

## IT IS RECOMMENDED IN THE FOLLOWING CASES:

- Infantile cerebral palsy – different forms.
- Meningomyelocele.
- Brain and spinal cord injuries.
- Muscle diseases connected with paresis and paralysis.
- Genetic disorders connected with paresis.
- Other diseases connected with paresis, paralysis and injuries regarding locomotor system.
- After traumatic brain injury, spine with spinal cord and lower limbs injury.



Stroller is equipped with Light Grey rims. However there is possibility to change the rim's colour to black.



## STANDARD EQUIPMENT:

Light aluminium frame | Fully foldable | Wheels with glass fiber, front PU swivel with direction | blockage, Rear inflatable with 2 separate wheel brace | Front shock absorbers | Ergonomic folding system | Seat unit with ClickClack™ system for easy plug in/out | 1 piece style height adjustable push handle | Reflective safety marks | Seat unit that can be used forward or rearward facing | Fully reclining backrest to make the stroller as a bed | Height adjustable backrest | Angle adjustable backrest | 5-points safety belts with height adjustment | Adjustable headrest | Anatomical lateral supports | Depth and width adjustable seat | Pneumatic tilt-in-space system TILT-IQ™ | Height and angle adjustable footrest | 5 points feet stabilizing belts with heel support and foot | abduction-adduction adjustments | Angle adjustable footplate | Soft upholstery with spine load reducing and ventilating | Soft front bar with reflective safety marks |



**ADVANTAGES:**



**Stable structure**

The product ensures adequate stability and safety of a Patient.



**Reflective elements**

The product has reflective elements to ensure your safety.



**Lightweight construction**

The product is made of lightweight and durable materials.

**24**

**24 months warranty**

We offer 24 months warranty on the product. We also provide after sales service.



**Easy to use**

Intuitive way of adjustment and easy assembly.




**Movable seat**

You can turn the seat easily so that the child may be facing towards caregiver or movement direction.



**Tilt-in-space**

It is possible to change the tilt angle of the seat.



Five-point straps for keeping safety.



Front swivel wheels for easy maneuverability.



The stylish and modern upholstery offers the highest comfort.

**COLORS**

**CITY UPHOLSTERY:**



City Grass  
Size 1



City Rose  
Size 2



City Coffee  
Size 3

**BREATHABLE UPHOLSTERY:**



Black  
Optional\*



Black & White  
Optional\*



Red  
Optional\*




Grass  
Optional\*



Sea  
Optional\*

\* Optional colors available on special request. Lead time ca. 45 days.  
Colors **BLACK/BLACK & WHITE/LATTE** are only available until the stocks last



Symbol	Measurement	Unit	Size 1	Size 2	Size 3
A	Seat width	[cm]	19 - 25	22 - 30	24 - 34
L	Seat depth	[cm]	21 - 26	27 - 35	34 - 42
C	Backrest height	[cm]	44 - 55	49 - 59	55 - 67
E	Backrest width	[cm]	23	27	32
F	Footrests length	[cm]	13 - 29	16 - 40	16 - 43
	Footer length	[cm]	20	20	20
	Push handle height	[cm]	100	100	100
X	Backrest tilt angle	[°]	90 - 180	90 - 180	90 - 180
Y	Seat tilt angle (forward-facing)	[°]	3 - 20	3 - 25	3 - 30
Y <sub>1</sub>	Seat tilt angle (rearward facing)	[°]	3 - 22	3 - 17	3 - 20
Z	Footplate tilt angle	[°]	90 - 180	90 - 180	90 - 180
	Max. user weight	[kg]	35	35	50
	Indicative user height*	[cm]	80 - 100	100 - 130	130 - 150
Product dimensions & weight		Unit	Size 1	Size 2	Size 3
	Width	[cm]	59	66	67
	Width after folding	[cm]	59	66	67
	Length	[cm]	90	90	90
	Length after folding	[cm]	90	90	92
	Height	[cm]	100	100	100
	Height after folding	[cm]	45	48	53
	Weight	[kg]	19,5	21	22,5

**Note:** Measurement tolerance range is +/- 1 cm / 1/2 in.  
 \*Height cannot be main factor for selection of product size



**URSUS HOME™** is the perfect stroller designed for daily rehabilitation of children with locomotor dysfunction.

Stable construction and an extensive system of supports and belts makes children feel very confidently and safely in **URSUS HOME™**.

Consequently, the small patient is able to take a meal, play, learn and perform the therapy exercises in seating or half-lying position. Caregivers can easily change the seat height. Angle adjustment of the backrest, seat, and footrest allows your child to rest in supine position.

A particular advantage that will be appreciated by therapist is infinitely adjustable tilt-in-space function.

Rich equipment and wide range of adjustment makes **URSUS HOME™** fit easily to the needs and dysfunctions of the child.

The compact construction of the stroller is perfect for small and tight spaces.

Thanks to the ClickClack™ system we can quickly detach the seat and attach it to the other products that are in our offer (**Ursus™**, **Ursus Home™**, **Ursus Active™**). The URSUS™ seat is available in 3 sizes.



## IT IS RECOMMENDED IN THE FOLLOWING CASES:

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.

## STANDARD EQUIPMENT:

Adjustable headrest | Anatomical lateral supports | 5-points safety belts with height adjustment | Height adjustable backrest | Angle adjustable backrest | Depth and width adjustable seat | Pneumatic tilt-in-space system tilt-iq™ | Height and angle adjustable footrest | 5 points feet stabilizing belts with heel support and foot | abduction-adduction adjustments | Soft upholstery with spine load reducing and ventilating | function | Soft front bar with reflective safety marks | Hi-low frame equipped with castors with brakes and push | handle



## ADVANTAGES:



### Stable structure

The product ensures adequate stability and safety of a Patient.



### Reflective elements

The product has reflective elements to ensure your safety.



### Lightweight construction

The product is made of lightweight and durable materials.



### 24 months warranty

We offer 24 months warranty on the product. We also provide after sales service.



### Easy to use

Intuitive way of adjustment and easy assembly.




### Movable seat

You can turn the seat easily so that the child may be facing towards caregiver or movement direction.



### Tilt-in-space

It is possible to change the tilt angle of the seat.



Five-point straps for keeping safety.



Easy height adjustment with a lever.



The stylish and modern upholstery offers the highest comfort.

## COLORS


### CITY UPHOLSTERY:



City Grass  
Size 1



City Rose  
Size 2



City Coffee  
Size 3

### BREATHABLE UPHOLSTERY:



Black  
Optional\*




Black & White  
Optional\*



Red  
Optional\*




Grass  
Optional\*



Sea  
Optional\*

\* Optional colors available on special request. Lead time ca. 45 days.  
Colors **BLACK/BLACK & WHITE/LATTE** are only available until the stocks last



Symbol	Measurement	Unit	Size 1	Size 2	Size 3
A	Seat width	[cm]	19 - 25	22 - 30	24 - 34
L	Seat depth	[cm]	21 - 26	27 - 35	34 - 42
C	Backrest height	[cm]	44 - 55	49 - 59	55 - 67
E	Backrest width	[cm]	23	27	32
F	Footrests length	[cm]	13 - 29	16 - 40	16 - 43
	Footer length	[cm]	20	20	20
	Push handle height	[cm]	86 - 96	86 - 96	86 - 96
X	Backrest tilt angle	[°]	90 - 180	90 - 180	90 - 180
Y	Seat tilt angle	[°]	0 - 25	0 - 25	0 - 25
Z	Footplate tilt angle	[°]	90-180	90-180	90-180
	Max. user weight	[kg]	35	35	50
	Indicative user height*	[cm]	80 - 100	100 - 130	130 - 150
Product dimensions & weight			Size 1	Size 2	Size 3
	Width	[cm]	55	55	55
	Length	[cm]	78	78	78
	Height (max)	[cm]	105	110	115
	Weight	[kg]	17	18,5	19,5

**Note:** Measurement tolerance range is +/- 1 cm / ½ in.  
 \*Height cannot be main factor for selection of product size.



**URSUS ACTIVE™** is a special wheelchair designed for daily rehabilitation of children with locomotor dysfunctions. Furthermore, this wheelchair may be very useful concerning orthopaedic disorders and lower limbs relief needed after surgery.

The wheelchair is equipped with anti-tippers wheels that prevent the patients to tip over backwards. Optional\*ly, **Ursus Active™** can also be equipped with push handles for a caregiver

Tilt-in-space function allows you to smoothly adjust the tilt angle of the seat which helps you to relief hips by moving the weight of body from pelvis area to the back.

Thanks to the ClickClack™ system we can quickly detach the seat and attach it to other products that are in our offer (**Ursus™, Ursus Home™, Ursus Active™, Quark™**).



## IT IS RECOMMENDED IN THE FOLLOWING CASES:

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.

## STANDARD EQUIPMENT:

Light aluminium frame | Large, inflatable rear wheels in sizes from 22" to 26" | equipped with quick-release couplings and strings | 2 independent rear brakes | Front shock absorbers | Seat unit with system ClickClack™ system for easy stick in/out | Fully reclining backrest to make the stroller as a bed (couch function) | Adjustable backrest | Anatomical lateral supports | 5-points safety belts with height adjustment | Height adjustable backrest | Angle adjustable backrest | Reflective safety marks | Pneumatic tilt-in-space system TILT-IQ™ | Height and angle adjustable footrest | 5 points feet stabilizing belts with heel support and foot | abduction-adduction adjustments | Angle adjustable footplate | Breathable upholstery | Soft front bar with reflective safety marks



## ADVANTAGES:



### Stable structure

The product ensures adequate stability and safety of a Patient.



### Reflective elements

The product has reflective elements to ensure your safety.



### Lightweight construction

The product is made of lightweight and durable materials.

# 24

### 24 months warranty

We offer 24 months warranty on the product. We also provide after sales service.



### Easy to use

Intuitive way of adjustment and easy assembly.



### Movable seat

You can turn the seat easily so that the child may be facing towards caregiver or movement direction.



### Tilt-in-space

It is possible to change the tilt angle of the seat.



Ergonomic push handle with height adjustment



Anti-tipping supports prevent the truck from tipping over.



Small front casters helps to maneuver even in small rooms.

## COLORS


### CITY UPHOLSTERY:



City Grass  
Size 1



City Rose  
Size 2



City Coffee  
Size 3

### BREATHABLE UPHOLSTERY:



Black  
Optional\*



Black & White  
Optional\*



Red  
Optional\*




Grass  
Optional\*



Sea  
Optional\*

\* Optional colors available on special request. Lead time ca. 45 days.  
Colors **BLACK/BLACK & WHITE/LATTE** are only available until the stocks last



Symbol	Measurement	Unit	Size 1	Size 2	Size 3
A	Seat width	[cm]	19 - 25	22 - 30	24 - 34
L	Seat depth	[cm]	21 - 26	27 - 35	34 - 42
C	Backrest height (without headrest)	[cm]	44 - 55	49 - 59	55 - 67
E	Backrest width	[cm]	23	27	32
F	Footrests length	[cm]	13 - 29	16 - 40	16 - 43
	Footer length	[cm]	20	20	20
	Push handle height	[cm]	85	85	85
X	Backrest tilt angle	[°]	90 - 180	90 - 180	90 - 180
Y	Seat tilt angle	[°]	0 - 25	0 - 25	0 - 25
Z	Footplate tilt angle	[°]	0 - 90	0 - 90	0 - 90
	Max. user weight	[kg]	35	35	50
	Indicative user height*	[cm]	80 - 100	100 - 130	130 - 150
Product dimensions & weight			Size 1	Size 2	Size 3
	Width	[cm]	62	66	70
	Length	[cm]	100	100	100
	Height	[cm]	87	92	87
	Weight	[kg]	18,5	20	22,5

**Note:** Measurement tolerance range is +/- 1 cm / 1/2 in.  
 \*Height cannot be main factor for selection of product size.



## Safety



**USS/USA\_001 LED lights**  
Front and rear LED lights for increased safety.



**USS\_002 Belt for moving the patients AM-P**  
Robust and ergonomic design as well as high durability facilitates daily care of the patient.



**USS/USA/USH\_003 Soft cover for tray**  
Soft cover protects the user from hitting the tray. Indicated for patients with impaired limbs coordination.



**USS\_004 Low head supporting panels**  
They provide additional head protection from unwanted contact with stroller's frame. Mounted around the push handle.



**USS\_022 Patient lift belt**  
Robust and ergonomic design as well as high durability facilitates daily care of the patient.



## Positioning



**USS/USA/USH\_101 Thigh abduction belts**  
They can be easily attached to the chair by use of Velcro straps. Their use allows to keep legs in abduction and also prevents the patient from slipping off of the device.



**USS/USA/USH\_107 Pelvic belt**  
Wide belt equipped with "HOLD AND PULL" system. It provides better stabilization in the device.



**USS/USA/USH\_113 Trunk stabilizing belt**  
It acts as an additional trunk protection which helps to keep a proper sitting position. The belt is fastened by use of Velcro straps and its size can be adjusted to the user.



**USS/USA/USH\_114 H Harness**  
It acts as a trunk stabilizer. HOLD&PULL™ system makes it easier to adjust the vest to the patient body. You can quickly detach the vest if necessary.



**USS/USA/USH\_121 Head-neck support**  
Perfectly stabilizes the head. Indicated for people with limited head control. You can mount it at different height.



**USS/USA/USH\_124 Loins support**  
It allows better stabilization of the lumbar back.



**USS/USA/USH\_125 6 points safety vest**  
Wide range adjustable vest that acts as a trunk stabilizer. HOLD&PULL™ system makes it easier to adjust the vest to the patient body.



**USS/USA/USH\_126 Trunk belt**  
Universal solution applicable to many of our products. It secures the user from falling out of the device. Anatomical shape and Velcro straps allow for a comfortable fit.

**U** URSUS™

**H** URSUS HOME™

**A** URSUS ACTIVE™

Akces-MED sp. z o.o.  
Jasionka 955B  
36-002 JASIONKA

export@akces-med.com  
+48 (17) 864 04 77  
www.akces-med.com



U  
A  
H

**USS/USA/USH\_127 Chest support**  
Width adjustable supports that provide additional side stabilization of the trunk. You can mount it at different height.



U  
A  
H

**USS/USA/USH\_128 Abduction block**  
Shaped abduction block provides with correct seating position. Can be easily removed.



U  
A  
H

**USS/USA/USH\_130 4 points safety belts**  
It acts as a trunk stabilizer. Wide range of adjustment helps to fit it to the user's needs. Equipped with HOLD&PULL™ system.



U  
A  
H

**USS/USA/USH\_131 Firm headrest made of foam**  
Made of poliurethane foam. Provides with basic support of the head



U  
A  
H

**USS/USA/USH\_134 Seat narrowing cushions 6 cm**  
Made of elastic foam trimmed with soft material. They are attached to the seat upholstery from both sides. Thanks to them you can narrow the seat by 6cm.



U  
A  
H

**USS/USA/USH\_137 Seat narrowing cushions 10cm**  
Made of elastic foam trimmed with soft material. They are attached to the seat upholstery from both sides. Thanks to them you can narrow the seat by 10cm.



U  
A  
H

**USS/USA/USH\_161 Sandals stabilizing feet and ankles**  
They protect and stabilize the ankle joint. They are made of elastic material, recommended for users with increased muscle tone. They can be used with or without shoes.



U  
A  
H

**USS/USA/USH\_162 Dynamic ankle stabilizers (FP-22)**  
Dynamic FlexPoint® ankle stabilizers provide great support when stabilizing the user's feet. Incorrect positioning of the feet, despite appropriate stabilization of the pelvis and chest, makes it impossible to keep the correct sitting position.



U  
A  
H

**USS/USA/USH\_163 Anatomic thigh abduction belts (FP-22)**



## BodyMap Positioning



U

**USS\_309 U cut seat cushion BodyMap® A-**  
Vacuum seat cushion allows for stabilisation and proper positioning of hips.



U

**USS\_310/314 U cut cotton cover for BodyMap®**  
Made of soft material it offers protection and absorbs moisture.



U

**USS\_311/315 U cut VisMemo cover for BodyMap®**  
Made of thermo-active foam for extra bed sore prevention. Allows for longer and more comfortable use of the cushion.



**USS\_312/316 U cut aerating-massaging pad for BodyMap®**  
Aerating-massaging pad to strengthen the prophylactic anti-bed sore function in patients using BodyMap. It can only be used together with the terry cover. The effect of being 'suspended in the air' creates a skin-friendly microclimate.

U URSUS™

H URSUS HOME™

A URSUS ACTIVE™

Akces-MED sp. z o.o.  
Jasionka 955B  
36-002 JASIONKA

export@akces-med.com  
+48 (17) 864 04 77  
www.akces-med.com



**USS\_313 Back cushion BodyMap® B**  
The vacuum back cushion provides comfort and correction for patients with scoliosis and poor posture. Cushion maintains and controls required posture in the stroller.



**USS\_343/344 Elastic cover for BodyMap®**



**USS\_350/351 3D mesh cover for BodyMap®**



### Comfort



**USS\_402 Umbrella**  
Attachable to the stroller's frame, perfectly protects from the sunlight or rain.



**USS/USA/USH\_403 Tray**  
Made of plastic. It supports the therapy and allows you to play, eat, and work.



**USS\_405 Folding canopy without side panels**  
Protects against both sun and rain.



**USS\_408 Rain Cover**  
An excellent protection from the rain, snow, and wind either for the user and stroller.



**USS\_409 Mosquito Net**  
Excellent protection against insects



**USS/USA/USH\_410 Cotton headrest cover**  
It helps to keep the headrest clean. Easy to assemble/disassemble, it can be washed.



**USS/USA/USH\_410 Cup holder**  
An universal holder mounted to the frame. It enables easy and quick access to a bottle of drink.



**USS/USA/USH\_414 Plexiglas® Tray**  
Made of Plexiglas®. It supports the therapy and allows you to play, eat, and work.



**USS\_415 Wheel covers (4 pcs.)**  
Use it on the dirt tires when you want to take the products indoor or transport in the trunk. .



**USS\_416 Summer footmuff**  
Protects legs and trunk. This esthetic body sheet fits the height of the user covering also the footrest of the pushchair.



**USS\_417 Winter footmuff**  
Protects legs and the torso area. A very tasteful cover of the body, adjusting to the height of the user, covering the footrest of the chair. Insulated lining provides extra protection against the cold.



**USS\_418 Hood for winter footmuff**  
Can be attached to the sleeping bag. Protects head against frost, wind and moisture.

**U** URSUS™

**H** URSUS HOME™

**A** URSUS ACTIVE™

Akces-MED sp. z o.o.  
Jasionka 955B  
36-002 JASIONKA

export@akces-med.com  
+48 (17) 864 04 77  
www.akces-med.com



**USS\_420 Folding canopy with side panels**  
Side panels protect not only against sun and rain but also the wind.



**USA\_426 Decorated spokes guard**  
Spokes guards 20", made OF polycarbonate 1,5mm, outsider diameter 415mm, insi de diameter 100mm



**USS/USA/USH\_429 Soft footrest cover**  
Comfortable pad, protecting the feet from direct contact with the footrest



**USS\_442 Front wheel covers**  
Use it on the dirt tires when you want to take the products indoor or transport in the trunk. .



**USS/USA\_447 Winter muffs**  
Comfortable and warm muffs protect your hands against a chill during walks.



**USS/USH/USA\_448 Lap tray**



## Storage



**USS/USA/USH\_503 Lady's bag**  
This capacious bag designed for mothers will be perfect for carrying all the necessities.



**USS/USA/USH\_504 Men's bag**  
Capacious men-style bag for carrying all the necessary things.



**USS\_505 Under seat storage basket**  
Enables you to carry small objects while you stroll. It is attached to the base of the stroller.



**USS\_506 Travel Bag Standard**  
Perfect solution to protect folded stroller against damage (e.g. during transport). Easily fits the shape and the size of a stroller.



**USS/USA/USH\_507 Oxygen tank case**  
Oxygen tank case can be attached to the stroller by use of Velcro straps with adjustable length. Made to measure to perfectly fit the patient's oxygen tank.



**USS/USA\_509 Organizer**  
A handy bag attached to the push handle. Specially designed compartments and pockets make it easy to keep order and have all the necessary items at hand. The organizer bag has two useful holes for a cup or a bottle of water.

**U** URSUS™

**H** URSUS HOME™

**A** URSUS ACTIVE™

Akces-MED sp. z o.o.  
Jasionka 955B  
36-002 JASIONKA

export@akces-med.com  
+48 (17) 864 04 77  
www.akces-med.com



## Build-up



U  
A

### USS/USA\_607 Hi-Low frame Ursus Home

Is a comfortable and functional solution in home conditions. Adjustable seat height greatly facilitates contact with peers, and push handles make maneuvering, even in tight spaces, much easier.



A  
H

### USA/USH\_608 Outdoor frame

Means tours and travels in very comfortable conditions either for a child and caregivers. It provides a convenient ride in varied terrain. Swivel front wheels ensure smooth changing of direction.



U  
H

### USS/USH\_609 Active frame

is a mode that gives your child a freedom of movement. Thoughtful solution enables independent movement of the user or with assistance of a caregiver. Now, therapy, learning, and fun will be even more enjoyable and effective.




U

### USS\_611 Kid-sit Buggypod PERLE



U

### USS\_612 Buggypod PERLE Adapter



U  
A  
H

### USA/USH\_615 URSUS seat



U  
A  
H

### USS/USA/USH\_623 QUARK frame

Will convert **URSUS™** seat into a comfortable rehabilitation chair. Innovative system of height adjustment ensures smooth and easy operation. Included push handle allows easy relocation of the device.