

ACTIVALL™

Dynamic standing frame



Class I Medical Device

**This is a medical device.
Use it in accordance with the user manual or the label.**

Description

As far this dynamic standing frame is concerned, it is highly beneficial to the process of daily rehabilitation which takes place not only in hospitals and rehabilitation wards but also in user's house. This device enables the user to stand and walk. Moreover, children are able to take first steps and adults can either keep or regain fitness.

Both light and stable frame is helpful in independent relocation in this device with slight involvement of the user. The correct adjustment of braces makes the user feel secure as well as comfortable during performing activities. Furthermore, this product may be used by users in different age groups due to a vast array of sizes.

ACTIVALL™ plays a significant role because of the fact that user has the possibility of gaining new experience and, more to the point, being independent as far as possible.

Indications

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use

Intuitive way of adjustment and easy assembly.



Independence

User can upright by himself independently.



2 in 1

The product has 2 utility functions for use in 1 single device. Dynamic and static function.



24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.

Colors



Size 1



Size 2



Size 3



Size 4



Size 5



Static function

Standing frame ACTIVALL™ with stabilizer of the lower limbs.



Dynamic function

Standing frame ACTIVALL™ with separator of the lower limbs.



Activall CROSS™

ACTIVALL CROSS™ makes it possible for the user to adopt upright position and move within house, hospital ward or outside.

Standard equipment

Adjustable chest support

Adjustable pelvic support

Adjustable hand grip






Pelvic harness


Separator of the lower limbs OR
Stabilizer of the lower limbs

Frame with castors equipped with brakes





Important measurements

Sizes	 Pelvic/chest support width	 Pelvic support height	 Chest support height	 Pelvic/chest support depth*	 Pelvic/chest support perimeter
Size 1	26 cm	44 - 60 cm	52 - 70 cm	22 cm	87 cm
Size 2	28 cm	66 - 82 cm	77 - 95 cm	25 cm	97 cm
Size 3	31 cm	77 - 102 cm	89 - 116 cm	27 cm	106 cm
Size 4	37 cm	83 - 112 cm	92 - 126 cm	29 cm	124 cm
Size 5	42 cm	98 - 120 cm	108 - 143 cm	35 cm	144 cm



Scan the QR code to see the full table.

*Measured without harness stabilizing the hips.
 **Height cannot be main factor for selection of product size.
 Note: Measurement tolerance range is +/- 1 cm / ½ in.
 Note: Approximate values, depending on the product configuration.

Sizes	 Indicative user height**	 Max. user weight
Size 1	70 - 100 cm	25 kg
Size 2	100 - 130 cm	40 kg
Size 3	110 - 150 cm	60 kg
Size 4	110 - 170 cm	80 kg
Size 5	140 - 200 cm	95 kg

Weight and dimensions

Size	Width	Length	Height	Weight
Size 1	59 cm	73 cm	54 cm	13,5 kg
Size 2	60 cm	84 cm	79 cm	16,5 kg
Size 3	63 cm	84 cm	90 cm	18 kg
Size 4	66 cm	91 cm	96 cm	20,5 kg
Size 5	74 cm	98 cm	112 cm	23 kg

Accessories



Belt for moving patients AM-P

Robust and ergonomic design as well as high durability facilitates daily care of the patient.



TANGO Wheel with fixed direction of movement

Noiseless, easy to maintain, does not leave marks, does not corrode.



Forearm supports

Perfect for people with limited control of the upper body. They provide an additional side stabilization.




Hip supporting harness

Attached to the pelvic support and acts as a pelvic girdle stabilizer which secures the user from falling out of the device.



Castor bumpers

Bumpers mounted to the stander wheels, protect the device from hitting obstacles.



Scan the QR code,
to see all the accessories.