

AKVOSEGO™

Bath chairs



Class I Medical Device

**This is a medical device.
Use it in accordance with the user manual or the label.**

Description

AKVOSEGO™ bath chair is simple and innovative solution that is helpful while taking care of disabled people. The bath chair consists of solid frame that is made of aluminium and, moreover, it protects from adverse effects of water. Furthermore, frame is equipped with foam covers that protect bathtub or kids bath tub from scratching. Unique construction of frame makes it possible to adjust the tilt angle of footrest and backrest.

Furthermore, it is also possible to change distance as well as seat's tilt angle towards floor. It is highly beneficial due to the fact that product is correctly adjusted to user's needs which is very helpful while e.g. washing little user's hair.

Concerning upholstery, it is made of impermeable net. **AKVOSEGO™** bath chair can be used together with **BodyMap®** cushion either in a room or during bath because all **BodyMap®** cushions are totally waterproof! The use of **BodyMap®** cushion makes the stabilization of patient more efficient during bath.

Its upholstery and frame are available in different colours – it depends on its size.

Indications

- cerebral palsy (CP)
- paralysis or paresis
- poliomyelitis (Heine-Medin disease)
- limb injuries or limb amputation
- limb deformities
- joint arthritis and injuries



Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use

Intuitive way of adjustment and easy assembly.



Waterproof

This construction is neutral for harmful water influence.



Lightweight construction

The product is made of light metal alloys which positively affects the low weight of the device.

Colors



Blue
Size 0




Green
Size 1



Orange
Size 2



Blue
Size 3



DRY FLOOR TECH™

DRY FLOOR TECH™. In all of our bath products the frame is filled with polyurethane foam which prevents water from getting inside.



Innovative articulated system

An innovative articulated system with a wide range of frame adjustments makes it possible to adjust the product to the user's measurements and needs.











Adjustable seat depth

An adjustable seat depth helps to fully adapt the product to the user's height, which provides a comfortable and stable sitting position.

Standard equipment



Important measurements

Sizes	 Product width	 Seat depth	 Headrest height	 Backrest height	 Seat height (from the floor)	 Lower leg length	 Indicative user height*	 Max user weight
Size 0	37 cm	23 - 28 cm	16 cm	31 cm	11 - 44 cm	26 cm	< 100 cm	25 kg
Size 1	42 cm	28 - 33 cm	16 cm	40 cm	11 - 44 cm	31 cm	100 - 120 cm	30 kg
Size 2	45 cm	32 - 40 cm	19 cm	44 cm	13 - 49 cm	36 cm	120 - 140 cm	50 kg
Size 3	50 cm	37 - 50 cm	23 cm	50 cm	15 - 54 cm	41 cm	130 - 170 cm	75 kg

Weight and dimensions

Size	Width	Width after folding	Length	Length after folding	Height	Height after folding	Weight
Size 0	42 cm	42 cm	70 cm	64 cm	80 cm	28 cm	4 kg
Size 1	42 cm	42 cm	75 cm	67 cm	85 cm	28 cm	4,5 kg
Size 2	45 cm	45 cm	91 cm	76 cm	96 cm	29 cm	5 kg
Size 3	50 cm	50 cm	97 cm	88 cm	113 cm	31 cm	8 kg

Scan the QR code to see the full table.




*Height cannot be the main parameter while choosing the product size.
Note: Measurement tolerance range is +/- 1 cm / ½ in.
Note: Approximate values, depending on the product configuration.

Accessories



Fixing kit for BodyMap® K cushion

A set of Velcro straps to fasten BodyMap cushions.



Thigh wedge

Shaped abduction block provides with correct seating position. Can be easily removed.



AKVO

Active waterproof platform for bath devices for disabled children.



BodyMap® K cushion

BodyMap® K is a stabilizing mattresses which can be easily shaped to the patient's body. Due to that quality, the BodyMap® products support every part of the patient's body and therefore are ideal in the bedsore prophylaxis.



Scan the QR code,
to see all the accessories.