

CAT I™

Static standing frame



MD Class I Medical Device

**This is a medical device.
Use it in accordance with the user manual or the label.**

Description

Stander **CAT I™** was created with a view to help children in adoption of the upright position. The upright standing contributes to the acquisition by children a new experience from the nearest environment and stimulates their psycho-physical development. In addition, it regulates the physiological processes and oxygenates the body of the child. There are many advantages during use of the **CAT I™**.

Adjustable tilt angle (from 25° to 85°) even with a child in the device. Extensive system of belts and supports makes the child feel safe and secure.

As far as this product is concerned, it is made of ecologic and environmentally friendly materials, and is also hand painted.

Indications

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use

Intuitive way of adjustment and easy assembly.



Eco-friendly materials

We use only human and environment friendly materials.



Pneumatic support

Gas spring helps in proper tilt angle adjustment of the device.



Quiet operation

Using of gas spring in the product ensures quiet operation.



3D foot buckles

Footrests with 3D adjustment allow for smooth foot setting in all directions.



Spacer upholstery

The device is equipped with breathable, hypoallergenic 3D spacer fabric. Now you can be sure that the skin stays dry for longer!




24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.

Colors

ATMO™ upholstery

Breathable / Antiallergic / Easy to clean / Reduced environment impact



Grass
Size 1



Bowl for toys

Tray with bowl for toys.



Gas spring

Gas spring helps you to set the proper tilt angle of the device.









3D foot adjustment



They allow for a smooth foot setting in all directions thanks to the use of the hinge ball.

Standard equipment



Important measurements

						
Size	Pelvic supports width	Chest supports width	Pelvic supports height adjustment	Chest supports height adjustment	Tilt angle	Pelvic/Chest supports depth
Size 1	12 - 23 cm	12 - 21 cm	38 - 60 cm	50 - 72 cm	25 - 85°	12 cm

		
Size	Indicative user height*	Max. user weight
Size 1	70 - 100 cm	25 kg

Weight and dimensions

Size	Width	Length	Height	Weight
Size 1	54 cm	100 cm	73 cm	17 kg

*Height cannot be the main parameter while choosing the product size.
Note: Measurement tolerance range is +/- 1 cm / 1/2 in.
Note: Approximate value, depending of product configuration.

Scan the QR code to see the full table.



Accessories



3D foot adjustment

They allow for a smooth foot setting in all directions thanks to the use of the hinge ball.



Head supporting belt

Attachable to the headrest by use of Velcro straps. Indicated for people with muscle weakness to prevent falling of the head.



Tray with bowl

It allows you to store even more toys or other things needed for therapy.



Steel wheel with brake

Steel wheel with brake, recommended when using wheels with a blockade of reverse movement on the stander or wheels with fixed direction of movement.



Scan the QR code,
to see all the accessories.