OMBRELO™

Special needs stroller

me



This is a medical device.

Use it in accordance with the user manual or the label.

Description

Light and easy to fold special needs stroller for children and adults. It has innovative seat **BodyPrint™** that was equipped with unique belts system. When adjusted perfectly gives the optimal positioning and comfort of use even for all day seating. Adjustable seat depth and backrest height for optimal comfort and correct position while seating. 30° fixed tilt promotes upper body positioning for those with limited trunk stability or upper body strength.

OMBRELO™ stroller is very comfortable and stable with excellent maneuverability.

Adjustable swing-away multi footplates for easy come out of the stroller FlipFlap $^{\text{m}}$ allow to adjust each foot separately. They are equipped with foot straps for 2 ways of stabilization.

The stroller is available in 4 sizes.

Size 5 is dedicated for adults with weight up to 120 kg.

OMBRELO™ passed the Crash Test in accordance with EN 7176/19.

Indications

- Cerebral palsy various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis SM).
- · Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- · Genetic disorders.
- Degeneration and neurological diseases.
- · Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use

Intuitive way of adjustment and easy assembly.



High durability

The construction ensures high durability of the device - maximum user weight is 120 kg!



24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.



Lightweight construction

The product is made of light metal alloys which positively affects the low weight of the device.



FlipFlap™ footrest

Special footrests system ensures a wide range of adjustment. Reclining platforms make getting on and off a device easier.



Colors

COMFORT™ upholstery



Moon Size 2



Loft Size 4



Sand Size 3



Onyx Size 5



Aloe

Optional for each size



Rose

Optional for each size

Waterproof upholstery - **HYDRO™**



Sea

Optional for each size



Sky

Optional for each size

COMFORT™ available while stocks last



Green grass Size 2



Brown Sugar Size 4



Hot Red Size 3



Grey Vulcano Size 5



Flip Flap™ footrest

Special footrests system ensures a wide range of adjustment. Reclining platforms make getting on and off a device easier.



BodyPrint™ upholstery

BodyPrint™ upholstery was equipped with unique belts system for seat and backrest. When adjusted perfectly gives the optimal positioning and comfort of use.




TransFix™ anchors

OMBRELO™ had been successfully crash tested and fulfilled all the requirements of standard EN ISO 7176/19. It makes possible to use the stroller in specially adjusted motor vehicles, equipped with the anchoring points and other needed equipment for forward-facing transport.

Standard equipment



Important measurements



Weight and dimensions

Size	Width	Width after folding	Length	Length after folding	Height	Height after folding	Weight
Size 2	57 cm	34 cm	125 cm	130 cm	101 cm	37 cm	18,5 kg
Size 3	63 cm	35 cm	137 cm	133 cm	104 cm	42 cm	18,5 kg
Size 4	67 cm	35 cm	137 cm	150 cm	104 cm	36 cm	20 kg
Size 5	72 cm	34 cm	137 cm	128 cm	108 cm	35 cm	20,5 kg



Scan the QR code to see the full table.

*Height cannot be the main parameter while choosing the product size.

Note: Measurement tolerance range is +/- 1 cm / ½ in.

Note: Approximate values, depending on the product configuration.

Accessories



Loins support

It allows better stabilization of the lumbar back.



Oxygen tank holder

Oxygen tank case can be attached to the stroller by use of Velcro straps with adjustable length. Made to measure to perfectly fit the patient's oxygen tank.



Umbrella

Attachable to the stroller's frame, perfectly protects from the sunlight or rain.



Trunk stabilizing belt

It acts as an additional trunk protection which helps to keep a proper sitting position. The belt is fastened by use of Velcro straps and its size can be adjusted to the user.



Safety rail

Makes the user safe in the stroller. Protects against falling out.



Scan the QR code, to see all the accessories.