

PARAPION™

Static standing frame



MD Class I Medical Device

**This is a medical device.
Use it in accordance with the user manual or the label.**

Description

We know very well how important it is for the user to adopt the standing position, that is why we have created a product called **PARAPION™** which is ideally suited for this purpose. **PARAPION™** allows you for full and safe upright standing.

Stable construction, ergonomic knee supports, wide pelvic and chest belts are the basic elements ensuring safety and comfort during the therapy.

The user, who is placed in this device, is able to perform rehabilitation treatment as well as simple domestic or professional activities.

The standing frame is equipped with four noiseless castors which enable easy relocation of the user.

Indications

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use


Intuitive way of adjustment and easy assembly.

24⁺


24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.


Colors



Size 1



Size 2



Size 3



Chest support

Ergonomic designed with wide cushions for excellent patient stabilization.
(available as accessory)



Independently adjustable knee supports

Adjustable and ergonomically shaped knee supports perfectly stabilize the lower limbs. (available as accessory)









Forearm supports

Perfect for people with limited control of the upper body. They provide additional side stabilization. (available as accessory)

Standard equipment





Important measurements

Sizes	 Pelvic/chest support width	 Pelvic/chest support circumference	 Chest support height	 Pelvic support height	 Knee supports width	 Knee supports height
Size 1	23 cm	56 - 76 cm	73 - 111 cm	61 - 90 cm	12 - 22 cm	18 - 45 cm
Size 2	31 cm	80 - 102 cm	92 - 126 cm	80 - 104 cm	16 - 24 cm	24 - 49 cm
Size 3	38 cm	95 - 124 cm	107 - 163 cm	92 - 122 cm	20 - 31 cm	31 - 71 cm



Scan the QR code to see the full table.

Size	 Indicative user height*	 Max user weight
Size 1	90 - 140 cm	30 kg
Size 2	110 - 170 cm	55 kg
Size 3	140 - 200 cm	95 kg

*Height cannot be the main parameter while choosing the product size.
 Note: Measurement tolerance range is +/- 1 cm / ½ in.
 Note: Approximate values, depending on the product configuration.

Weight and dimensions

Size	Width	Length	Height	Weight
Size 1	56 cm	75 cm	73 cm	17 kg
Size 2	56 cm	87 cm	92 cm	21,5 kg
Size 3	67 cm	93 cm	107 cm	27,5 kg

Accessories



Belt for moving patients AM-P

Robust and ergonomic design as well as high durability facilitates daily care of the patient.



Pelvic support from Paramobil

Ergonomic designed with wide cushions for excellent patient stabilization.



Chest support from Paramobil

Ergonomic designed with wide cushions for excellent patient stabilization.



TANGO Wheel with brake

Noiseless, easy to maintain, does not leave marks, does not corrode.



Tray

Made of wood. It supports the therapy and allows you to play, eat, and work.



Scan the QR code,
to see all the accessories.