

URSUS™ & URSUS HOME™

Special needs strollers



URSUS™



URSUS HOME™

This is a medical device.
Use it in accordance with the user manual or the label.

Description

URSUS™

Ultralight stroller **URSUS™** is designed for daily rehabilitation of children with locomotor disabilities. As far as its construction is concerned, the frame and seat are made of high-strength aluminum, and the chassis is additionally cushioned. **URSUS™** is equipped with **ClickClick™** feature, which means that you can easily swap the seats and adjust the stroller to the most suitable mode.

The device has a wide range of adjustments such as width, seat depth, backrest height, backrest tilt angle, infinitely adjustable tilt-in-space function – **TILT-IQ™**, reversible seat. **URSUS™** has also been equipped with spatial, breathable upholstery that prevents excessive sweating as well as bedsores, and more to the point, positively influences the comfort of use. Moreover, the stroller can easily be converted into a lounger so that the child can sleep in it.

The product is available in 3 different sizes with a very wide range of accessories that will help you to meet the individual needs of user.

Indications

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use

Intuitive way of adjustment and easy assembly.



Lightweight construction

The product is made of light metal alloys which positively affects the low weight of the device.



Reflective elements

The product has reflective elements to ensure your safety.



Tilt-In-Space

The device is able to change hip angle.



Movable seat

You can turn the seat easily so that the child may be facing towards caregiver or movement direction.

24⁺

24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.

Reversible
seat






Tilt angle
adjustment

Colors

Waterproof upholstery **HYDRO™**



CADET
Standard equipment

REVERSIBLE COTTON-WASHABLE PAD COMFORT




CADET
Standard equipment



RUBY*
Standard equipment




OCEAN
Standard equipment



SPACE
Standard equipment



LAVA
Standard equipment



LILAC
Standard equipment

*Available until the stocks last.




Safety

5-point belts provide the highest level of security.



Upholstery

Reversible cotton-washable pad COMFORT.












Footrest

Height and angle adjustable footrest

Standard equipment



Important measurements

Sizes	 Seat width	 Seat depth	 Backrest height	 Backrest width	 Backrest tilt angle	 Footplate tilt angle	 Footplate length	 Indicative user height*	 Max user weight
Size 2	22 - 30 cm	27 - 35 cm	49 - 59 cm	27 cm	90 - 180°	90 - 180°	16 - 40 cm	100 - 130 cm	35 kg
Size 3	24 - 34 cm	34 - 42 cm	55 - 67 cm	32 cm	90 - 180°	90 - 180°	16 - 43 cm	130 - 150 cm	50 kg

Weight and dimensions

Size	Width	Width after folding	Length	Length after folding	Height	Height after folding	Weight
Size 2	63 cm	63 cm	90 cm	90 cm	100 cm	48 cm	21 kg
Size 3	67 cm	67 cm	90 cm	90 cm	100 cm	53 cm	22,5 kg



Scan the QR code to see the full table.

*Height cannot be the main parameter while choosing the product size.

In extreme positions the handle can collide with the backrest when changing its angle. In such cases, the handle should be lifted while it is being adjusted.

Note: Measurement tolerance range is +/- 1 cm / ½ in.

Note: Approximate values, depending on the product configuration.

Description

URSUS HOME™

URSUS HOME™ is the perfect stroller designed for daily rehabilitation of children with locomotor dysfunction.

Stable construction and an extensive system of supports and belts makes children feel very confidently and safely in **URSUS HOME™**.

Consequently, the small patient is able to take a meal, play, learn and perform the therapy exercises in seating or half-lying position. Caregivers can easily change the seat height. Angle adjustment of the backrest, seat, and footrest allows your child to rest in supine position.

A particular advantage that will be appreciated by therapist is infinitely adjustable tilt-in-space function.

Rich equipment and wide range of adjustment makes **URSUS HOME™** fit easily to the needs and dysfunctions of the child.

The compact construction of the stroller is perfect for small and tight spaces.

Thanks to the ClickClack™ system we can quickly detach the seat and attach it to the other products that are in our offer (**URSUS™**, **URSUS HOME™**).

The **URSUS™** seat is available in 3 sizes.

Indications

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use

Intuitive way of adjustment and easy assembly.



Lightweight construction

The product is made of light metal alloys which positively affects the low weight of the device.



Reflective elements

The product has reflective elements to ensure your safety.



Tilt-In-Space

The device is able to change hip angle.



Movable seat

You can turn the seat easily so that the child may be facing towards caregiver or movement direction.


24⁺

24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.

Colors

Waterproof upholstery **HYDRO™**



CADET
Standard equipment

REVERSIBLE COTTON-WASHABLE PAD COMFORT




CADET
Standard equipment



RUBY*
Standard equipment



OCEAN
Standard equipment



SPACE
Standard equipment



LAVA
Standard equipment



LILAC
Standard equipment

*Available until the stocks last.




Safety

5-point belts provide the highest level of security.



Upholstery

Reversible cotton-washable pad COMFORT.












Footrest

Height and angle adjustable footrest

Standard equipment



Important measurements


Sizes	 Seat width	 Seat depth	 Backrest height	 Backrest width	 Backrest tilt angle	 Seat tilt angle	 Footrests length	 Indicative user height*	 Max user weight
Size 2	22 - 30 cm	27 - 35 cm	49 - 59 cm	27 cm	90 - 180°	90 - 180°	16 - 40 cm	100 - 130 cm	35 kg
Size 3	24 - 34 cm	34 - 42 cm	55 - 67 cm	32 cm	90 - 180°	90 - 180°	16 - 43 cm	130 - 150 cm	50 kg

Weight and dimensions

Size	Width	Length	Height (maximal)	Weight
Size 2	55 cm	78 cm	110 cm	18,5 kg
Size 3	55 cm	78 cm	115 cm	19,5 kg

Weight and dimensions of the HI-LOW frame

Size	Width	Length	Height (without push handles)	Weight
Universal	55 cm	78 cm	47 - 58 cm	10,6 kg



Scan the QR code to see the full table.

*Height cannot be the main parameter while choosing the product size.
Note: Measurement tolerance range is +/- 1 cm / ½ in.
Note: Approximate values, depending on the product configuration.

Accessories **URSUS™** / **URSUS HOME™**



PA-VM-15 cushion

VISmemo™ memory foam used in the cushion increases the seating comfort. A bolster at the front of the cushion helps maintain the position of the pelvis, reducing the possibility of sliding off the seat.



Oxygen tank holder

Oxygen tank case can be attached to the stroller by use of Velcro straps with adjustable length. Made to measure to perfectly fit the patient's oxygen tank.




Umbrella

Attachable to the stroller's frame, perfectly protects from the sunlight or rain.



Travel Bag Standard

An excellent solution to protect the folded stroller from damages (e.g. during the transportation). It perfectly fits the stroller shape and dimensions.



Scan the QR code,
to see all the accessories.