

# JORDI™ & JORDI HOME™



**JORDI™**  
Positioning chair



**JORDI HOME™**  
Special needs stroller

**This is a medical device.  
Use it in accordance with the user manual or the label.**

## Description

### JORDI™

---

**JORDI™** is a school-therapeutic chair, which works great for therapy and for fun, as well as other daily activities such as learning or eating meals. At home, kindergarten, school or in a group of peers, **JORDI™** will be an irreplaceable help tool in the rehabilitation process, both in the aspects of movement and social integration.

Durable design and regulation of the tilt-in-space function are the main advantages of the seat, which is especially appreciated by therapists. Furthermore, the seat has height, width and depth adjustment. Already in the standard it becomes the ideal tool in the rehabilitation process and a long list of accessories will allow us to fit the product to individual needs.

Another advantage is the medical upholstery, easy to clean and disinfect. For more demanding customers we prepared a medical upholstery called **EXTREME-MED™**, which has certified non-flammability (CRIB5), is antibacterial and antifungal.

Available in **3** pretty, vibrant colors.

Made of safe materials for the user. Available in **4** sizes.

## Indications

---

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



## Special features



### **Stable structure**

The device ensures adequate stability and safety of user.



### **Easy to use**

Intuitive way of adjustment and easy assembly.



### **Eco-friendly materials**

We use only human and environment friendly materials.



### **Tilt-In-Space**

The device is able to change hip angle.



### **Hold & Pull**

A special system used in our vests and hip belts. The ergonomic D-shaped grip makes it easy for the user to adjust the product quickly and comfortably. The system is designed also for people with hand dysfunctions.



### **24-month warranty**

We offer 24-month warranty on the product. We also provide after sales service.

## Colors

Medical upholstery **MEDICUS™**




**MED SEA**  
Standard



**MED SKIN**  
Optional

Medical upholstery **EXTREME-MED™**



**PILLARBOX\***  
Optional

\*Available while stocks last



## Seat depth adjustment

Seat depth adjustment and backrest tilt angle available in standard equipment.



## Tilt-In-Space

---

Tilt is a space function is available thanks to the legs height adjustment.



## Medical upholstery

---






Medical upholstery used in this product is vapor permeable, waterproof and easy to disinfect.


## Standard equipment

---





Important measurements

					
Sizes	Seat width	Seat depth	Backrest height	Floor to seat height	Backrest tilt angle
Size 1	25 cm	23 - 29 cm	33 cm	26 - 30 cm	90 - 135°
Size 2	30 cm	26 - 34 cm	39 cm	31 - 39 cm	90 - 135°
Size 3	37 cm	30 - 40 cm	46 cm	37 - 49 cm	90 - 135°
Size 4	45 cm	35 - 45 cm	55 cm	42 - 56 cm	90 - 135°



Scan the QR code to see the full table.

		
Size	Indicative user height*	Max user weight
Size 1	90 - 110 cm	50 kg
Size 2	100 - 120 cm	65 kg
Size 3	110 - 150 cm	80 kg
Size 4	150 - 180 cm	120 kg

\*Height cannot be the main parameter while choosing the product size.  
Note: Measurement tolerance range is +/- 1 cm / ½ in.  
Note: Approximate values, depending on the product configuration.

Weight and dimensions

Size	Width	Length	Height	Weight
Size 1	37 cm	48 cm	56 cm	6,5 kg
Size 2	42 cm	54 cm	68 cm	8 kg
Size 3	52 cm	60 cm	80 cm	10 kg
Size 4	64 cm	68 cm	92 cm	12 kg

## Description

# JORDI HOME™

---

**JORDI HOME™** is a special needs stroller which is perfect both for therapy and recreation, as well as other every day activities such as learning or having meals. The stroller can be used at home, school, in the kindergarten or in a group of peers. **JORDI HOME™** turns out to be an indispensable tool in the rehabilitation process, both in terms of therapy and social integration.

A solid construction and the regulation of the tilt-in-space function are the main advantages of the stroller, which are especially appreciated by the therapists. The seat in the **JORDI HOME™** stroller has an adjustable height, width and depth. The HexagonPro™ adjustable headrest guarantees a correct and extremely comfortable head support. The 4-point safety vest additionally enhances the comfort and protection of the user.

The **JORDI HOME™** stroller is an ideal tool in the rehabilitation process with the standard equipment, however, the wide range of accessories allows the stroller to be fully adjusted to the individual needs of the user. The stroller is equipped with four swivel wheels with a brake, which is why it moves quietly and is exceptionally maneuverable. The user's feet are secured on the adjustable footrest.

Another advantage is a medical upholstery, easy to clean and disinfect. For more demanding customers we prepared a medical upholstery called **EXTREME-MED™**, which has certified non-flammability (**CRIB5**), is antibacterial and antifungal. The upholstery is available in three pretty, vibrant colours.

**JORDI HOME™** is made of user friendly materials. It is available in 4 sizes.

## Indications

---

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.





## Special features



### Stable structure

The device ensures adequate stability and safety of user.



### Easy to use

Intuitive way of adjustment and easy assembly.



### Eco-friendly materials

We use only human and environment friendly materials.



### Tilt-In-Space

The device is able to change hip angle.



### Hold & Pull

A special system used in our vests and hip belts. The ergonomic D-shaped grip makes it easy for the user to adjust the product quickly and comfortably. The system is designed also for people with hand dysfunctions.



### 24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.

## Colors

Medical upholstery **MEDICUS™**




**MED SEA**  
Standard



**MED SKIN**  
Optional

Medical upholstery **EXTREME-MED™**



**PILLARBOX\***  
Optional

\*Available while stocks last



## Seat depth adjustment

Seat depth adjustment and backrest tilt angle available in standard equipment.



## Tilt-In-Space

---

Tilt is a space function is available thanks to the legs height adjustment.



## Medical upholstery


---

Medical upholstery used in this product is vapor permeable, waterproof and easy to disinfect.


## Standard equipment





## Important measurements



Sizes	Seat width	Seat depth	Backrest height	Floor to seat height	Backrest tilt angle
Size 1	25 cm	23 - 29 cm	33 cm	32 - 36 cm	90 - 135°
Size 2	30 cm	26 - 34 cm	39 cm	39 - 47 cm	90 - 135°
Size 3	37 cm	30 - 40 cm	46 cm	54 - 57 cm	90 - 135°
Size 4	45 cm	35 - 45 cm	55 cm	50 - 64 cm	90 - 135°



Scan the QR code to see the full table.

Size	Indicative user height*	Max user weight
Size 1	90 - 110 cm	50 kg
Size 2	100 - 120 cm	65 kg
Size 3	110 - 150 cm	80 kg
Size 4	150 - 180 cm	120 kg

\*Height cannot be the main parameter while choosing the product size.  
 Note: Measurement tolerance range is +/- 1 cm / ½ in.  
 Note: Approximate values, depending on the product configuration.

## Weight and dimensions

Size	Width	Length	Height	Weight
Size 1	37 cm	77 cm	79 cm	9,5 kg
Size 2	42 cm	84 cm	94 cm	11,5 kg
Size 3	52 cm	100 cm	109 cm	13,5 kg
Size 4	64 cm	115 cm	124 cm	16 kg

# Accessories

---



## Hip guide

Gives great stabilization in hip area. Can narrow the seat width.

---



## Knee abductor

Adjustable and easy to remove. Prevent in legs crossing. Made of soft foam, upholstered with easy to wash medical fabric.


---



## Abduction block

Shaped abduction block provides with correct seating position. Can be easily removed.

---



## H vest

It acts as a trunk stabilizer. HOLD&PULL™ system makes it easier to adjust the vest to the patient body. You can quickly detach the vest if necessary.

---



## Laterall support

Width adjustable supports that provide additional side stabilization of the trunk. You can mount it at different height.

---

Scan the QR code,  
to see all the accessories.



JORDI™



JORDI HOME™