KIDOO[™] & KIDOO HOME[™] ●



This is a medical device. Use it in accordance with the user manual or the label.

Description **KIDOO**™

KIDOO™ is a therapeutic chair that is perfect for therapy and fun as well as other daily activities such as learning or taking meals. The device can be used at home, school, kindergarten or in a group of peers. **KIDOO™** turns out to be an indispensable tool in the rehabilitation process, both in terms of therapy and social integration.

This orthopedic chair provides the child with proper support and safety, moreover, it can be easily adjusted to meet your child's needs and dysfunctions. Pelvic belt height adjustable seat and armrests further enhance comfort and safety of the user.

It is made of user-friendly materials. Available in 5 sizes.

Indications

- Cerebral palsy various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



Special features



Stable structure The device ensures adequate stability and safety of user.



Eco-friendly materials We use only human and environment friendly materials.



24

Tilt-In-Space The device is able to change hip angle.

Intuitive way of adjustment and easy



Hold & Pull

A special system used in our vests and hip belts. The ergonomic D-shaped grip makes it easy for the user to adjust the product quickly and comfortably. The system is designed also for people with hand dysfunctions.

24-month warranty

Easy to use

assembly.


We offer 24-month warranty on the product. We also provide after sales service.

Colors

Medical upholstery **MEDICUS™**







Backrest tilt angle adjustment

Seat depth adjustment and backrest tilt angle available in standard equipment.



Tilt-In-Space



Medical upholstery

Medical upholstery used in this product is vapor permeable, waterproof and easy to disinfect.

Tilt is a space function is available thanks to the legs height adjustment – click for more information

Standard equipment





23 cm

26,5 cm

31 cm

35 cm

42 cm

• •


Size 1

Size 2

Size 3

Size 4

Size 5



26 cm

31 cm

35 cm

40 cm

44 cm

19,5 - 26 cm

24 - 31,5 cm

27,5 - 36,5 cm

33 - 44 cm

36 - 49 cm

90 - 116°

90 - 116°

90 - 116°

90 - 116°

90 - 116°

18 - 23,5 cm

22,5 - 30 cm

24,5 - 34,5 cm

31 - 41 cm

36 - 47 cm



Scan the QR code to see the full table.

			
Size	Indicative user height*	Max user weight	
Size 1	90 - 100 cm	20 kg	
Size 2	90 - 115 cm	30 kg	
Size 3	100 - 130 cm	40 kg	
Size 4	110 - 150 cm	55 kg	
Size 5	120 - 160 cm	70 kg	

*Height cannot be the main parameter while choosing the product size. Note: Measurement tolerance range is +/-1 cm $/\frac{1}{2}$ in. Note: Approximate values, depending on the product configuration.

Weight and dimensions

Size	Width	Length	Height	Weight
Size 1	34 cm	37,5 cm	48 cm	5 kg
Size 2	37 cm	41 cm	56 cm	6,5 kg
Size 3	42 cm	53 cm	63 cm	8 kg
Size 4	47 cm	58 cm	75 cm	11 kg
Size 5	55 cm	66 cm	81 cm	14 kg

Description **KIDOO HOME**™

KIDOO HOME™ is a special needs stroller which is perfect both for therapy and recreation, as well as other every day activities such as learning or having meals. The stroller can be used at home, school, in the kindergarten or in a group of peers. **KIDOO HOME™** turns out to be an indispensable tool in the rehabilitation process, both in terms of therapy and social integration.

Its construction is based on the therapeutic chair for children **KIDOO™**, and just like in the **KIDOO™** chair, the child is provided with a proper support and safety. Apart from that, the stroller can be easily adjusted to meet your child's needs and disfunctions.

The 4-point safety vest additionally enhances the comfort and protection of the user. The stroller has an adjustable seat and armrests height as well as the backrest tilt ankle. The **HexagonPro™** adjustable headrest guarantees a correct and extremely comfortable head support.

KIDOO HOME[™] stroller is set on a mobile base equipped with four swivel wheels with a brake, which is why it moves quietly and is exceptionally maneuverable. The user's feet are secured on the adjustable footrest with 2-point feet stabilizing belts.

KIDOO HOME™ is made of user friendly materials. It is available in 5 sizes.

Indications

- Cerebral palsy various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



Special features



Stable structure The device ensures adequate stability and safety of user.



Intuitive way of adjustment and easy assembly.



Eco-friendly materials We use only human and environment friendly materials.



24

Tilt-In-Space The device is able to change hip angle.



Hold & Pull

A special system used in our vests and hip belts. The ergonomic D-shaped grip makes it easy for the user to adjust the product quickly and comfortably. The system is designed also for people with hand dysfunctions.


24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.

Colors


Medical upholstery **MEDICUS™**





Backrest tilt angle adjustment

Seat depth adjustment and backrest tilt angle available in standard equipment.



Tilt-In-Space

Tilt is a space function is available thanks to the legs height adjustment – click for more information



Medical upholstery

Medical upholstery used in this product is vapor permeable, waterproof and easy to disinfect.


Standard equipment



Important measurements

42 cm

Size 5



44 cm

36 - 49 cm

90 - 116°



Scan the QR code to see the full table.

Size	Indicative user height*	Max user weight	
Size 1	90 - 100 cm	20 kg	
Size 2	90 - 115 cm	30 kg	
Size 3	100 - 130 cm	40 kg	
Size 4	110 - 150 cm	55 kg	
Size 5	120 - 160 cm	70 kg	

36 - 47 cm

*Height cannot be the main parameter while choosing the product size. Note: Measurement tolerance range is +/- 1 cm / $\frac{1}{2}$ in. Note: Approximate values, depending on the product configuration.

Weight and dimensions

Size	Width	Length	Height	Weight
Size 1	34 cm	37,5 cm	48 cm	5 kg
Size 2	37 cm	41 cm	56 cm	6,5 kg
Size 3	42 cm	53 cm	63 cm	8 kg
Size 4	47 cm	58 cm	75 cm	11 kg
Size 5	55 cm	66 cm	81 cm	14 kg

Accessories



Foot and ankle stabilizing sandals

They protect and stabilize the ankle joint. They are made of elastic material, recommended for users with increased muscle tone. They can be used with or without shoes.



Knee abductor

Adjustable and easy to remove. Prevent in legs crossing. Made of soft foam, upholstered with easy to wash medical fabric.



Narrowing pads 6/10 cm

Made of elastic foam trimmed with medical upholstery. They are attached to the armrests from both sides. Thanks to them you can narrow the seat by 6/10cm.



Laterall support

Width adjustable supports that provide additional side stabilization of the trunk. You can mount it at different height.

Scan the QR code, to see all the accessories.



