

MAMALU™

Special needs stroller



Class I Medical Device

**This is a medical device.
Use it in accordance with the user manual or the label.**

Description

Modern and practical with unique design – **MAMALU™** is a stroller created for the youngest users, which will work perfectly in all conditions. Its construction is based on a light and durable umbrella frame which means that it folds easily and quickly, and its small dimensions make it easy to transport after folding.

MAMALU™ grows together with your child due to the wide range of adjustments the stroller is equipped with. The headrest and the 5-point safety belts can be mounted at different heights, the seat has an adjustable depth and there are three levels of height adjustment of the footrest. Using the stroller is great pleasure. Independent shock absorbers and swivel PU wheels make the **MAMALU™** stroller extremely maneuverable. Comfort and convenience are provided by a reversible pad **COMFORT**, which is covered with soft, skin-friendly **MINKY** fabric on one side, and **FLEX** waterproof fabric on the other.

Check the **MAMALU™** stroller during a walk and see that it is created for you and your child.

Indications

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use

Intuitive way of adjustment and easy assembly.



Lightweight construction

The product is made of light metal alloys which positively affects the low weight of the device.



Reflective elements

The product has reflective elements to ensure your safety.


24⁺

24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.

Colors

FLEX upholstery



CADET
Standard equipment

Reversible cotton-washable pad comfort




CADET
Optional




RUBY*
Optional




OCEAN
Optional



SPACE
Optional



LAVA
Optional



LILAC
Optional

*Available until the stocks last.



The umbrella frame

The umbrella frame of the stroller makes it fold easily and quickly, and its small dimensions make it easy to transport after folding.



Grow with your child

Adjustable headrest, 5-point safety belts, footrest and seat depth make the stroller grow with your child.



Reversible pad COMFORT









Reversible pad COMFORT – a combination of soft, skin-friendly MINKY fabric with waterproof upholstery.



Standard equipment




Important measurements

Sizes	 Seat width*	 Seat depth	 Backrest width*	 Backrest height*	 Backrest tilt angle (relative to the seat)	 Footrest length
Size 1	38 (32) cm	22, 29 cm	38 (32) cm	57 cm	95°	12, 20, 28, 29 cm
Size	 Indicative user height**	 Max user weight				
Size 1	90 - 130 cm	50 kg				

Weight and dimensions

Size	Width	Width after folding	Length	Length after folding	Height	Height after folding	Weight
Size 1	61 cm	32 cm	87 - 94 cm	120 cm	111 cm	30 cm	12 kg

Scan the QR code
to see the full table.



*High rigid trunk supports narrow the seat / backrest by 6 cm to 32 cm.
 **Height cannot be the main parameter while choosing the product size.

Note: Measurement tolerance range is +/- 1 cm / ½ in.
 Note: Approximate values, depending on the product configuration.

Accessories



High lateral supports

They provide additional side support of the trunk. Also, they protect the user against the contact with the stroller frame.



Footmuff

Protects legs and the torso area. A very tasteful cover of the body, adjusting to the height of the user.



PA-VM-15 cushion

VISmemo™ memory foam used in the cushion increases the seating comfort. A bolster at the front of the cushion helps maintain the position of the pelvis, reducing the possibility of sliding off the seat.




Folding canopy

Protects against both sun and rain.



Reversible cotton-washable pad COMFORT

Reversible pad COMFORT – a combination of soft, skin-friendly MINKY fabric with waterproof upholstery.



Scan the QR code, to see all the accessories.