

# TATALU™

Special needs stroller



Class I Medical Device

**This is a medical device.  
Use it in accordance with the user manual or the label.**

## Description

---

**TATALU™** – a stroller which makes walks with your child even more comfortable. A new version of the stroller with a lightweight and durable construction designed to perform well in all conditions. Its umbrella frame folds easily and quickly and takes up little space after folding. A number of improvements enhance the comfort of use both for the child and caregiver. The stroller can be fully customized – you decide which accessories will work for you. **TATALU™** – means comfort and stabilization! A number of improvements, such as the **SOFT** upholstery covering the sides of the stroller to provide protection of the user, adjustable 5-point safety belts, headrest or footrest with the possibility of changing the height, make the young user feel comfortable and safe. The **SOFT** upholstery meets the requirements of caregivers, which is confirmed by the **CIGARETTE TEST** proving that it is made of a non-flammable fabric. The increased comfort is due to an independent, reinforced shock absorption system and swivel front wheels which provide great maneuverability of the stroller. **TATALU™** – see for yourself that walks with your child can be even more enjoyable!

## Indications

---

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



## Special features



### **Stable structure**

The device ensures adequate stability and safety of user.



### **Easy to use**

Intuitive way of adjustment and easy assembly.



### **Lightweight construction**

The product is made of light metal alloys which positively affects the low weight of the device.



### **Reflective elements**

The product has reflective elements to ensure your safety.

24<sup>+</sup>

### **24-month warranty**

We offer 24-month warranty on the product. We also provide after sales service.


## Colors

---

Upholstery **SOFT**



**GALAXY**



**GRAFEN**



**CORAL**



**GOLD**



**STONE**



**UNIVERSE**



## Umbrella frame

---

The umbrella frame of the stroller makes it fold easily and quickly, and its small dimensions make it easy to transport after folding.



## SOFT upholstery

---

Comfortable SOFT upholstery is made of smooth, skin-friendly material, which makes travelling in the stroller very enjoyable.



## Safety

---







5-point safety belts with height adjustment and side supports guarantee safety and proper stabilization while using the stroller.





## Standard equipment



## Important measurements

Sizes	 Seat width*	 Seat depth*	 Backrest width	 Backrest height	 Backrest tilt angle (relative to the seat)	 Footrest length
Size 1	38 (30) cm	22, 30 cm	38 (30) cm	60 cm	95°	12, 20, 28, 39 cm

Size	 Indicative user height**	 Max user weight
Size 1	90 - 130 cm	50 kg

## Weight and dimensions

Size	Width	Width after folding	Length	Length after folding	Height	Height after folding	Weight
Size 1	62 cm	36 cm	104 cm	120 cm	104 cm	28 cm	13 kg

Scan the QR code  
to see the full table.



\*High rigid trunk supports narrow the seat / backrest by 8 cm to 30 cm.  
\*\*Height cannot be the main parameter while choosing the product size.

Note: Measurement tolerance range is +/- 1 cm / ½ in.  
Note: Approximate values, depending on the product configuration.



# Accessories



## Adjustable headrest

Adjustable height and proper shape of the headrest guarantee a perfect adjustment to the user's needs.



## Abduction block

They provide additional side support of the trunk. Also, they protect the user against the contact with the stroller frame.



## Folding canopy

Protects against both sun and rain.




## Dynamic ankle stabilizers (FP-22)

Dynamic FlexPoint® ankle stabilizers provide great support when stabilizing the user's feet.



## Tray

ABS thermoformable tray, perfect for recreation and during rehabilitation. It is also ideal for work and having meals. Tray is also available in light grey color.



Scan the QR code, to see all the accessories.