

JUMBO™

Positioning chair



Class I Medical Device

**This is a medical device.
Use it in accordance with the user manual or the label.**

Description

It can be safely used by children with increased muscle tone, as well as very flaccid children with decreased muscle tone. It is perfect for a lot of daily activities such as rehabilitation exercises, eating and playing, among others.

The medical upholstery **MEDICUS™** used in the chair is vapour-permeable, waterproof and easy to disinfect. Thus, it can be successfully used at home, individually, or in a hospital ward or rehabilitation center where it will be used by many people.

The foldable footrest is very helpful when putting your child on the seat. It can be disassembled for easy storage and transportation. The universal design of the footrest plate allows it to be mounted in two settings, which increases the range of footrest length adjustments.

Out of concern for people and the environment, the product is made of eco-friendly, safe materials.

The **JUMBO™** chair is especially recommended for people who have difficulty maintaining a sitting position on their own (especially with neurological and orthopedic conditions)

Indications

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



Special features



Stable structure

The device ensures adequate stability and safety of user.



Eco-friendly materials

We use only human and environment friendly materials.

24⁺

24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.



Easy to use

Intuitive way of adjustment and easy assembly.



Medical upholstery

Made of medical upholstery which is vapour-permeable and waterproof. It's safe for the users and easy to disinfect. It doesn't rub off, which makes it perfect to be used by different people, e.g. in rehabilitation centers. Ideal for use in areas of high risk of infection.



Hold & Pull

A special system used in our vests and hip belts. The ergonomic D-shaped grip makes it easy for the user to adjust the product quickly and comfortably. The system is designed also for people with hand dysfunctions.

Colors

Medical upholstery **MEDICUS™**



MED SEA
Standard



Independently adjustable side supports

Independently adjustable side supports ensure proper trunk stabilization.



Available adjustments

Due to available adjustments we are able to set the desired seat depth and tilt angle of the backrest.









Foldable footrest

The foldable footrest with extended adjustment range provides adequate foot stabilization.

Standard equipment





Important measurements

						
Sizes	Seat width (with tray)	Seat depth	Backrest height	Backrest width	Backrest tilt angle	Footplate length
Size 1	38 (35) cm	18 - 30 cm	41 cm	38 cm	90 - 115°	20 cm
Size 2	44 (41) cm	19 - 38 cm	45 cm	44 cm	90 - 110°	24 cm
Size 3	49 (46) cm	29 - 55 cm	61 cm	49 cm	90 - 115°	27 cm



Scan the QR code to see the full table.

		
Size	Indicative user height*	Max user weight
Size 1	< 110 cm	30 kg
Size 2	100 - 150 cm	45 kg
Size 3	110 - 180 cm	75 kg

*Height cannot be the main parameter while choosing the product size.
 Note: Measurement tolerance range is +/- 1 cm / ½ in.
 Note: Approximate values, depending on the product configuration.

Weight and dimensions

Size	Width	Length (Length with folded footrest)	Height	Weight
Size 1	49 cm	84 (74) cm	84 cm	14 kg
Size 2	54 cm	103 (90) cm	95 cm	20 kg
Size 3	59 cm	118 (106) cm	105 cm	24,5 kg

Accessories



HexagonPro™ headrest

Individual "ladyfingers" on the Hexagon Pro™ are adjusted to cushion and support any head position. It has been designed to create the most comfort for user by supporting and protecting their head.



Knee abductor

Adjustable and easy to remove. Prevent in legs crossing.



VisoPro™ headrest

The individual parts of the headrest help protect the head from unwanted rotation and from falling to the sides. The upper part supports the head laterally, the middle part controls its position around the ear area, while the lower part supports the area around the base of the skull.



4 points safety vest

Acts as a trunk stabilizer. Wide range of adjustment helps to fit it to the user's needs. Equipped with HOLD&PULL™ system.



Tray

Tray made of ecological wood. It lets patient play and supports therapy. Ideal for work and meals.



Scan the QR code, to see all the accessories.