

ULISES EVO™

Special needs stroller



Class I Medical Device

**This is a medical device.
Use it in accordance with the user manual or the label.**

Description

ULISES EVO™ is a perfect stroller designed for daily rehabilitation of children and youth with locomotor dysfunctions. It is particularly useful for users requiring additional stabilization.

The stroller has adjustable seat depth and backrest tilt angle which greatly facilitates its functionality.

With simple and light construction it can be easily folded for transporting purposes. Furthermore, relatively small size allows you to fit **ULISES EVO™** into an undersized car trunk.

Stable structure of the stroller, belts and supports system as well as reflective elements make the user feel safe. The backrest can be easily bend back and adjustable push handle increases the stroller maneuverability.

Rich equipment and wide range of adjustments make **ULISES EVO™** fit easily to the individual needs.

Indications

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.



Special features



Stable structure

The device ensures adequate stability and safety of user.



Reflective elements

The product has reflective elements to ensure your safety.



Lightweight construction

The product is made of light metal alloys which positively affects the low weight of the device.



Easy to use

Intuitive way of adjustment and easy assembly.



24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.

Colors

Waterproof upholstery - **HYDRO™**



Space gray



Neptune*

*Optional color. Lead time may be extended up to 30 days.



Wide range of adjustments

Wide range of adjustments make ULISES EVO™ fit easily to the individual needs.



High rigid trunk supports

High rigid trunk supports.



Footrest












Footrest with height and tilt angle adjustments, 2-point belts fastening the feet.



Standard equipment



Important measurements

Sizes	 Seat width without stabilizing side supports	 Seat width with stabilizing side supports	 Width between footrest angle adjusters	 Seat depth	 Headrest height	 Backrest width without stabilizing side supports	 Backrest width with stabilizing side supports	 Backrest tilt angle	 Footrests tilt angle	 Indicative user height*	 Max user weight
Size 1	31 cm	21 cm	27 cm	22 - 29 cm	40 - 70 cm	28 cm	21 cm	90 - 150°	90 - 180°	80 - 120 cm	55 kg
Size 1a	31 cm	21 cm	27 cm	29 - 37 cm	43 - 80 cm	28 cm	21 cm	90 - 150°	90 - 180°	100 - 150 cm	55 kg
Size 2	37 cm	27 cm	33 cm	29 - 37 cm	43 - 80 cm	34 cm	27 cm	90 - 150°	90 - 180°	100 - 150 cm	75 kg
Size 2a	37 cm	27 cm	33 cm	35 - 45 cm	51 - 90 cm	34 cm	27 cm	90 - 150°	90 - 180°	130 - 180 cm	75 kg
Size 3	43 cm	32 cm	39 cm	35 - 45 cm	51 - 90 cm	40 cm	32 cm	90 - 150°	90 - 180°	130 - 180 cm	75 kg

Weight and dimensions

Size	Width	Width after folding	Length min.	Length max.	Length medium	Length after folding	Height	Height after folding	Weight
Size 1	58 cm	58 cm	82 cm	130 cm	94 cm	104 cm	102 cm	40 cm	18,5 kg
Size 1a	58 cm	58 cm	98 cm	145 cm	110 cm	107 cm	102 cm	45 cm	19,5 kg
Size 2	64 cm	64 cm	98 cm	145 cm	110 cm	107 cm	102 cm	45 cm	20 kg
Size 2a	64 cm	64 cm	102 cm	165 cm	127 cm	114 cm	102 cm	50 cm	21,5 kg
Size 3	70 cm	70 cm	102 cm	165 cm	127 cm	114 cm	102 cm	50 cm	22 kg



Scan the QR code to see the full table.

*Height cannot be the main parameter while choosing the product size.

Note: In extreme positions of the handle, it may collide with the backrest when changing its angle. In such cases, the handle should be raised while adjusting.

Note: Measurement tolerance range is +/- 1 cm / ½ in.

Note: Approximate values, depending on the product configuration.

Accessories



Trunk stabilizing belt

It acts as an additional trunk protection which helps to keep a proper sitting position.



Loins support

It allows better stabilization of the lumbar back. Can be mounted at an appropriate height.



Abduction block

It is recommended for people with increased adductor tension as it separates lower limbs and helps maintain correct pelvic alignment while sitting.



Stroller cover 'n carry bag

Protects against both sun and rain.



Scan the QR code,
to see all the accessories.