

# RACER BODYMAP®

Special needs stroller



RACER BODYMAP®



RACER EVO BODYMAP®

**This is a medical device.  
Use it in accordance with the user manual or the label.**

## Description

# RACER BodyMap®

---

**RACER BodyMap®** is a perfect solution for people with lack or limited body stabilization.

**RACER BodyMap®** is a stroller that enables carrying children and adults with limited mobility. It can be used either in the user's home or rehabilitation centres. Stable construction as well as system of vacuum cushions BodyMap® make the user feel safe and comfortable.

**BodyMap®** system consists of:

- **BodyMap® A+** seat cushion with pelvis stabilization and abduction block – its profiled shape contributes to the correct sitting position and the abduction block prevents crossing of the legs.
- **BodyMap® B+** backrest cushion with trunk side stabilization – provides additional trunk stabilization perfectly adapting to the shape of the spine. Used both for posture alignment and comfort increase. You must no longer use extra "fillers" and other rollers.
- **BodyMap® D headrest** – it makes your head appropriately supported. It easily adapts to the shape of the head. The headrest is filled with granulate whose properties allow additional stabilization from every side of the head.
- **Lower limbs separator** with calf support belt helps to stabilize and keep the lower limbs in correct position, additionally prevents from crossing of the legs.
- **Electric pump** – greatly facilitate shaping of the cushions. Most importantly, you are able to shape the cushions without necessity of removing the user from the stroller. In order to change its shape again, you should inflate some air inside.

## Indications

---

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.





## Special features



### **Stable structure**

The device ensures adequate stability and safety of user.



### **Reflective elements**

The product has reflective elements to ensure your safety.



### **BodyMap® system**

Stroller is designed to fully stabilize the user using the vacuum BodyMap®.



### **Easy to use**

Intuitive way of adjustment and easy assembly.



### **24-month warranty**

We offer 24-month warranty on the product. We also provide after sales service.

## Colors



Black



## BodyMap® system

A unique advantage of RACER BodyMap® stroller is connection between its construction and BodyMap® vacuum system



## Calf belt stabilizer

---

The calf belt stabilizer perfectly stabilizes the user's legs and protects them from possible injury.



## Wheels

---

Rear wheels with shock absorbers provide comfort even on bumpy terrain.

## Standard equipment





## Description

# RACER EVO BodyMap®

---

**RACER EVO BodyMap®** is a perfect solution for people with lack or limited body stabilization.

**RACER EVO BodyMap®** is a stroller that enables carrying children and adults with limited mobility. It can be used either in the user's home or rehabilitation centres. Stable construction as well as system of vacuum cushions BodyMap® make the user feel safe and comfortable.

**BodyMap®** system consists of:

- **BodyMap® A+** seat cushion with pelvis stabilization and abduction block – its profiled shape contributes to the correct sitting position and the abduction block prevents crossing of the legs.
- **BodyMap® B+** backrest cushion with trunk side stabilization – provides additional trunk stabilization perfectly adapting to the shape of the spine. Used both for posture alignment and comfort increase. You must no longer use extra" fillers" and other rollers.
- **BodyMap® D headrest** – it makes your head appropriately supported. It easily adapts to the shape of the head. The headrest is filled with granulate whose properties allow additional stabilization from every side of the head.
- **Lower limbs separator** with calf support belt helps to stabilize and keep the lower limbs in correct position, additionally prevents from crossing of the legs.
- **Electric pump** – greatly facilitate shaping of the cushions. Most importantly, you are able to shape the cushions without necessity of removing the user from the stroller. In order to change its shape again, you should inflate some air inside.

## Indications

---

- Cerebral palsy – various types and forms.
- Brain and spinal cord injuries with lower limbs paralysis and weakness, balance and motor coordination disorders with different etiology.
- Demyelinating disease (e.g. multiple sclerosis – SM).
- Meningocele.
- Muscles diseases with weakness and paralysis (dystrophy, myasthenia gravis, myositis).
- Genetic disorders.
- Degeneration and neurological diseases.
- Other condition connected with muscles weakness, paralysis and musculoskeletal disorders.
- After head and spinal cord injuries and spine trauma with paralysis and paresis.





## Special features



### **Stable structure**

The device ensures adequate stability and safety of user.



### **Easy to use**

Intuitive way of adjustment and easy assembly.



### **Reflective elements**

The product has reflective elements to ensure your safety.



### **24<sup>+</sup> 24-month warranty**

We offer 24-month warranty on the product. We also provide after sales service.



### **BodyMap® system**

Stroller is designed to fully stabilize the user using the vacuum BodyMap®.



## Colors



Black



## BodyMap® system

A unique advantage of RACER BodyMap® stroller is connection between its construction and BodyMap® vacuum system



## Calf belt stabilizer

---

The calf belt stabilizer perfectly stabilizes the user's legs and protects them from possible injury.



## Wheels

---

Rear wheels with shock absorbers provide comfort even on bumpy terrain.

## Accessories **RACER BODYMAP® / RACER EVO BODYMAP®**



### Side covers

Side panels protect not only against sun and rain but also the wind.



### Tray

ABS thermoformable tray, perfect for recreation and during rehabilitation. It is also ideal for work and having meals.



### PA-VM-15 cushion

VISmemo™ memory foam used in the cushion increases the seating comfort. A bolster at the front of the cushion helps maintain the position of the pelvis, reducing the possibility of sliding off the seat.



### Folding canopy

Protects against both sun and rain.



### Mosquito net

Excellent protection against insects.



Scan the QR code,  
to see all the accessories.

## Important measurements

---

The biggest advantage of our cushions is their perfect fit to the user body. Every cushion is made individually according to the dimensions given in order form, which allows to meet every user's needs.

**Note:** In extreme positions of the handle, it may collide with the backrest when changing its angle. In such cases, the handle should be raised while adjusting.

Each **BodyMap**® product is made by hand with greatest care and attention to every detail, therefore, the waiting time for individual stroller is up to 30 days.

Ask for details in your shop or call the **Akces-MED** consultant **+48 516 183 883**.

To learn more about the **BodyMap**® system, please visit [www.bodymapsystem.com](http://www.bodymapsystem.com).

Scan the QR code  
to see the full table.

