

HIPPO+™ NEW VERSION

Special needs stroller



Class I Medical Device

**This is a medical device.
Use it in accordance with the user manual or the label.**

Description

We are presenting new and improved special stroller **HIPPO+™** which was equipped with many interesting and useful solutions. It is still perfect stroller for daily care with supreme positioning features.

Equipped in standard with swivel front wheels and reversible seat unit for your child to face you or the world. Bigger adjustable canopy as well as leg cover protect from wind, rain and sun and, more to the point, they enable rest irrespective of weather conditions. Shock absorbers soften the vibration during riding and 5 points safety belts connected with side supports provide your child with correct position and safety.

Reach standard and optional equipment help you customize the stroller for your child's needs. It is highly recommended for the youngest up to 2 years.

Indications

- Infantile cerebral palsy – different forms.
- Meningomyelocele.
- Brain and spinal cord injuries.
- Muscle diseases connected with paresis and paralysis.
- Genetic disorders connected with paresis.
- Other diseases connected with paresis, paralysis and injuries regarding locomotor system.
- After traumatic brain injury, spine with spinal cord and lower limbs injury.



Special features



Stable structure

The device ensures adequate stability and safety of user.



Easy to use

Intuitive way of adjustment and easy assembly.



Lightweight construction

The product is made of light metal alloys which positively affects the low weight of the device.



Reflective elements

The product has reflective elements to ensure your safety.

24⁺

24-month warranty

We offer 24-month warranty on the product. We also provide after sales service.



Movable seat

You can turn the seat easily so that the child may be facing towards caregiver or movement direction.

Colors



MELANGE*



VULCANO*



LOTUS



AQUA*



DENIM*



FLAMING*

*The lead time for selected colors (DENIM/FLAMING/AQUA/VULCANO/MELANGE) may be longer.



Movable seat

The seat can be put in the forward or rearward position.



Standard equipment

A proper stabilization of the user is provided by a wide range of standard equipment.



Standard equipment



+ Rain cover




+ Leg cover



+ Seat and backrest narrowing cushions





Important measurements



Sizes	Seat width	Seat depth	Backrest width	Backrest height	Backrest tilt angle	Footrests length	Footplate length	Footrests tilt angle
-------	------------	------------	----------------	-----------------	---------------------	------------------	------------------	----------------------

Size 1	20*/31 cm	24 cm	31 cm	40*/45 cm	90 - 150°	13 - 25 cm	14 cm	90 - 180°
--------	-----------	-------	-------	-----------	-----------	------------	-------	-----------

Size 2	34 cm	28 cm	34 cm	55 cm	90 - 150°	15 - 27 cm	18,5 cm	90 - 180°
--------	-------	-------	-------	-------	-----------	------------	---------	-----------

Size	Indicative user height**	Max user weight
------	--------------------------	-----------------

Size 1	70 - 100 cm	25 kg
--------	-------------	-------

Size 2	90 - 120 cm	35 kg
--------	-------------	-------

Weight and dimensions

Size	Width	Width after folding	Length	Length after folding	Height	Height after folding	Weight
Size 1	64 cm	64 cm	118 cm	112 cm	105 cm	56 cm	19 kg
Size 2	67 cm	67 cm	120 cm	112 cm	113 cm	56 cm	21 kg

Scan the QR code to see the full table.



*Measurement starts from narrowing pad.

**Height cannot be the main parameter while choosing the product size.

Note: Measurement tolerance range is +/- 1 cm / ½ in.
Note: Approximate values, depending on the product configuration.

Accessories



Loins support

It allows better stabilization of the lumbar back. Can be mounted at an appropriate height.



Reduction set consisting of additional seat

The cushion narrows down the seat by 6 cm. It is made of elastic foam trimmed with soft material.



Headrest with belt

Mounted in place of the standard headrest. It can be precisely adjusted to the user's needs. Soft head belt additionally stabilizes the forehead.



Foot and ankle stabilizing sandals

They protect and stabilize the ankle joint. They can be used with or without shoes.



Tray

ABS thermoformable tray, perfect for recreation and during rehabilitation. It is also ideal for work and having meals.



Scan the QR code, to see all the accessories.